

V.Mueller®

Cardiovascular Instruments


Instruments cardiovasculaires

Instrumente für die Kardiovaskularchirurgie

Strumenti chirurgici cardiovascolari

Instrumentos quirurgicos cardiovasculares

Instrumentos cirurgicos cardiovascular



V. Mueller Cardiovascular Instruments

Introduction	4
Elevators.....	14
Rib Shears.....	16
Retractors	19
Dissectors	31
Scissors	32
Needle Holders	46
Forceps.....	58
Clamps	76
Miscellaneous / Suction Tubes	125
Index	132

Dear customers,

More than 100 years ago in Chicago, Vinzenz Mueller started a small storefront business to provide physicians and surgeons with fine instruments and exceptional customer service. Mr. Mueller was guided by the principle that keeping all promises made to his customers - even if at a cost - was essential to the long-term success of his business. Since then, V. Mueller® has set the standard for service in the industry, which is still part of our company culture today.

The tradition of making fine surgical instruments continued when William Merz started as an instrument apprentice for Mr. Mueller in 1928. Driven by the idea of improving the lives of patients, Mr. Merz used his specialized knowledge and skills to translate surgeon problems or ideas into more efficient surgical instruments. His passion for advancing quality surgical instrumentation helped build his reputation as one of the best-known instrument makers in the world.

We are thankful for the thousands of innovative surgeons, operating room and sterile processing professionals, and other healthcare personnel who have collaborated with us to develop our high-quality, innovative new products. Today, the V. Mueller and Snowden-Pencer™ instrument lines include more than 25,000 instruments used in hospitals across the globe.

We are also proud of our craftsmen and supporting staff who continue our tradition of providing high-quality surgical instrumentation and customer service for customers everywhere.



William Merz



Reputation

Our standards of quality originated with our founder, Vinzenz Mueller. In 1893, he brought the knowledge of German instrument craftsmanship to the United States and established a tradition that has carried through to today. His protégé, William Merz, established the V. Mueller name as one synonymous with innovation and collaboration.

In 1951, Leonard Snowden and George Pencer introduced the revolutionary tungsten carbide insert to the jaws of surgical instruments. CareFusion is now proud to offer the Snowden-Pencer family of products to our customers.

These men dedicated their lives to improving surgical instruments and patient care. Our customers have always expected high quality and first-rate technology from our surgical instruments, which are promises CareFusion continues to deliver on today.

Commitment

At CareFusion, our commitment to quality stems from our ability to understand customer needs and exceed industry expectations.

Our sales force has more than 125 knowledgeable V. Mueller representatives, some with more than 30 years of experience. Also, our V. Mueller and Snowden-Pencer customer support team ensures accurate order processing and technical support. We are all here to help.

Once our products are in your hands, our service quality does not stop. We offer many opportunities for continuing education. Plus, all of our surgical instruments are backed by a lifetime warranty against defects in materials and workmanship.

William Merz started the V. Mueller tradition of working closely with well-recognized surgeons to develop surgical instruments used all over the world.

These surgeons have included Adson, Atkinson, Bailey, Balfour, Barraquer, Baumgartner, Beckman, Castaneda, Cooley, Cosgrove, Cottle, Crile, Cushing, Davis, Deakey, Derf, Finochietto, Glassman, Hibbs, Jackson, Jacobsen, Jannetta, Kelly, Lalonde, Mars, Mayo, McCulloch, McPherson, Merz, Metzenbaum, Mueller, Paparella, Pierse, Potts, Ramirez, Rhoton, Rochester, Spetzler, Tebbetts and Webster.



Quality

CareFusion prides itself on the quality craftsmanship of our V. Mueller and Snowden-Pencer surgical instruments. Though technology has progressed over the years, all of our instruments maintain a handcrafted component.

We know that small details make a real difference in the hands of surgeons. That is why we consider our surgical instruments an art form and attention to detail critical. This mindset and pride ensure our surgical instruments are among the highest in quality.

Our entire manufacturing process reflects our dedicated specializations through:

- Training instrument craftsmen by apprenticeships
- Inspecting instruments for quality throughout manufacturing

- Placing instruments under several stages of heat treatment, followed by quick cooling to achieve the proper hardness and strength
- Deep cleaning instruments, and preventing their corrosion by passivation

We have also proudly maintained a manufacturing presence in the United States for more than 115 years. With state-of-the-art manufacturing and operational facilities in St. Louis, MO; Tucker, GA; Cleveland, OH; and Tuttlingen, Germany, our quality instruments are manufactured under consistent standards by highly skilled individuals.



“We give the best to our customers by doing a microstructure evaluation of the raw materials. This helps ensure that what our customers are getting in their hands is a good product.”

Joe Torok
Principal Scientist/Metallurgist
CareFusion

Innovation

The V. Mueller and Snowden-Pencer brands have a history of developing recognized innovations, including:

- The tungsten carbide insert process, enabling a superior grasping technique
- Snowden-Pencer laparoscopic hand instruments, incorporating the finest instrument jaw patterns with ergonomically contoured handles to better fit surgeon hands

CareFusion continues to embody this philosophy of innovation. We deeply understand our customers' needs by working closely with our surgical and clinical advisory boards.

Our research and development team of engineers uses modern computer-aided design, analysis tools and specialized test laboratories to produce high-quality, innovative designs. They then hand over these validated designs to our manufacturing team. These experienced craftsmen incorporate handcrafting with advanced machinery such as CNC mills and lathes to consistently deliver some of the highest-quality instruments on the market. In effect, every new product we release is built with the same superior standards that our customers have come to expect for the last 115+ years.



V. Mueller product warranty

We warranty each medical device and accessory bearing the V. Mueller and Snowden-Pencer brand names to be free of functional defects in workmanship and materials when used normally for its intended use. Any V. Mueller product proven defective will be replaced or repaired at no charge. Repairs or modifications performed by unauthorized personnel may void all product warranties, or affect device performance and efficacy.

CareFusion employees are our own toughest critics. From reviewing production at every stage to demanding test procedures in our development and production processes, we insist on providing quality instruments you can trust. That is why all of our instruments bearing the V. Mueller name carry a lifetime warranty and provide personalized support from our V. Mueller sales representatives and customer support specialists.

The Snowden-Pencer Products Guarantee

All Snowden-Pencer instruments are guaranteed to be free of defects in material and workmanship. In addition, the tungsten carbide insert is guaranteed against wear under normal use for five years on Diamond-Jaw, Micro Diamond-Jaw, Diamond-Points and Diamond-Edges instruments.

Diamond-SerEdge scissors are guaranteed against wear for three years.

Diamond-Edges Supercut and Diamond-Tip instruments are guaranteed against wear for one year.

After this guarantee period, service and maintenance are available for a nominal charge.



We have persevered for over a century to build and maintain a reputation for quality - and are proud to stand behind our work with a lifetime warranty. In 2001, the Mayo Clinic reported a broken instrument to their V. Mueller representative. As the representative inspected the instrument, he found an interesting mark: a date-stamp from October 12, 1955. The instrument was nearly half a century old. But bearing the V. Mueller name, it was covered by a lifetime warranty against defects in materials and workmanship. So, we replaced it at no charge - just another example of the V. Mueller reputation for quality and commitment.

Return goods policy

In general, returns are NOT authorized for:

- Any item older than 90 calendar days from the invoice date
- Any item purchased on a "special order" or custom basis
- Used, processed or sterilized instruments, apparatuses or equipment
- Expired or non-expired items with missing or damaged labeling or packaging
- Sterile merchandise without all manufacturer seals intact
- Vibratory-etched products
- Color-coded merchandise

Customer support

Authorized representatives

V. Mueller products are sold only by authorized V. Mueller representatives. All instruments and equipment of our manufacture bear the V. Mueller trademark.

For more information, please contact your sales representative, send an email to GMB-VM-Europe@carefusion.com or visit carefusion.com/vmueller.



How to use this catalog

Use this key to easily understand the information throughout this catalog, as indicated in the sample pages below:

- A.** This catalog is divided by surgical specialty and instrument type.
- B.** Each page has its own ruler marked in inches, centimeters and millimeters.
- C.** Category names appear in English, French, German, Italian, Japanese, Portuguese and Spanish.
- D.** A diamond symbol ♦ denotes the most popular instrument within a category.
- E.** An oval (CH2461) around a catalog number denotes that item is pictured.
- F.** **titanium** denotes titanium instruments.
- G.** Illustrations of tip and jaw patterns appear next to the instrument.

Notes

- The alphabetical and numerical indexes at the end of this catalog can help you locate instruments by name, type or catalog number.
- All instruments are surgical stainless steel unless otherwise specified.

Needle Holders		Forceps	
<p>Vital® WANGENSTEEN Needle Holder Straight jaws. Tungsten Carbide inserts. 3,600 teeth per square inch. Overall length 10-3/4in (27.3cm).</p> <p>CH2433</p>		<p>Delicate Touch® Micro Forceps Overall length 8-1/4in (21cm).</p> <p>CH8030-21 Straight tip.</p> <p>CH8035-21 Curved tip.</p>	
<p>Vital® FINOCHIETTO Needle Holders Tungsten Carbide inserts. 3,600 teeth per square inch.</p> <p>CH2461 Overall length 8in (20.3cm).</p> <p>CH2471 Overall length 10-1/2in (26.7cm).</p>		<p>Delicate Touch® Micro Tying Ring Tip Forceps Straight jaws, length 17.5cm</p> <p>CH8016-17 Overall length 6-3/4in (17.5cm).</p> <p>CH8018-21 Overall length 8-1/4in (21cm).</p>	
<p>Vital® Microsurgery Needle Holders Spring handle. Smooth Tungsten Carbide inserts.</p> <p>CH2452 Straight jaws. Flat lock. Overall length 5-3/4in (14.6cm).</p> <p>CH2448 Curved jaws. Flat lock. Overall length 5-5/8in (14.3cm).</p> <p>CH2454 Straight jaws. Without lock. Overall length 5-3/4in (14.6cm).</p> <p>CH2457 Straight jaws. Flat lock. Overall length 7-1/4in (18.4cm).</p>		<p>Vital® CASTROVIEJO Needle Holders Spring handle. Flat lock. Straight jaws. Tungsten Carbide inserts. 3,600 teeth per square inch.</p> <p>CH2456 Overall length 7-1/4in (18.4cm).</p> <p>CH8589 Overall length 8-3/4in (22.2cm).</p>	
<p>Delicate Touch® Bayonet-Shaped Forceps Straight.</p> <p>CH8034-19 Overall length 7-5/8in (19.4cm).</p>		<p>Titanium Delicate Touch® GERALD Forceps Straight. Flat Handle. Platform tips.</p> <p>CH03760 Overall length 7in (17.8cm).</p> <p>CH03762 Overall length 8in (20.3cm).</p> <p>titanium</p>	

P4-34 carefusion.com/vmueller V. Mueller® Products and Services 1.800.323.9088 P4-19

How to measure instruments

The following measuring definitions can help you choose the right instruments for your needs:


Length

Length is measured in a straight line from the bottom of the instrument to the tip farthest away (A). On curved and angled instruments, the length is determined by a line drawn from the bottom of the instrument to the tip farthest away (B).

Jaw length

In clamps, the usable portion of the jaw is measured in a straight line to the bottom of the serrations (C). On partial occlusion clamps, the jaw length and jaw depth are measured as illustrated (D).

The dimensions listed are accurate as of this catalog’s publication date. Slight variations may occur due to metal crafting and finishing operations. Jaw and tip illustrations are approximately actual size unless otherwise noted.



Conversion table	
Inches to centimeters	Multiply inches (as decimals) by 2.54
Centimeters to inches	Multiply centimeters by 0.3937
Millimeters to inches	Multiply millimeters by 0.03937


Inches	Millimeters	Inches	Millimeters
1/16	1.588	9/16	14.288
1/8	3.175	5/8	15.875
3/16	4.736	11/16	17.463
1/4	6.350	3/4	19.050
5/16	7.938	13/16	20.638
3/8	9.525	7/8	22.225
7/16	11.113	15/16	23.813
1/2	12.700	1	25.400




Standard French (Charriere) gauge

The standard French, or Charriere, scale (Fr) is generally used in the size calibration of catheters and other tubular instruments. Based on the metric system, each unit is approximately 0.33 mm, differing 0.33 mm in diameter between consecutive sizes.

Example: 15 Fr indicates a diameter of 5 mm, and 30 Fr indicates a diameter of 10 mm.



Instrument angles



Protect your surgical instrumentation with Genesis™ containers

Protect your valuable surgical instrument investment with the V. Mueller brand Genesis container system - the name you trust in the OR. The Genesis container system was designed with hard-coated anodized aluminum to protect and maintain the sterility of your surgical instrumentation. Available in a wide range of sizes, the Genesis container system offers a variety of baskets and accessories to help you organize and speed the set-up of your instruments. The Genesis container system consistently provides a return on your investment and a green alternative to sterilization wrap.



Forceps/Elevators


Forceps/Élévateurs
Zangen/Heber
Pinze/Elevatori
Fórceps/Elevadores
Fórceps/Elevadores
鉗子/エレベーター

Chest Forceps

Long jaw, length 6in (15.2cm).

CH9300

Straight.
Overall length
12-1/2in (31.7cm).



CH9302

Curved jaw.
Overall length
11-7/8in (30.2cm).

CH9304

C-Shape.
Overall length
10-7/16in (26.5cm).

CH9306


Curved-on-side.
Overall length
12in (30.5cm).

CH9308

Curved handle.
Overall length
12-1/8in (30.8cm).

CH9312

45° Ring jaw.
Overall length
12-3/8in (31.4cm).



Chest Forceps

Short jaw, length
4in (10.2cm).

CH9314

Straight.
Overall length
12-1/2in (31.7cm).

CH9316

Curved jaw.
Overall length
12-1/8in (30.8cm).

CH9318

C-shaped.
Overall length
10-3/4in (27.3cm).

CH9320

Curved-on-side.
Overall length
12in (30.5cm).

CH9322

Curved handle.
Overall length
11-7/8in (30.2cm).



DOYEN Costal Elevators

Overall length 7in (17.8cm).

CH35

Right.
Blade 1in (2.5cm).

CH36

Left.
Blade 1in (2.5cm).

CH30

Right.
Blade 1-1/4in (3.2cm).

CH31

Left.
Blade 1-1/4in (3.2cm).



MATSON Rib Strippers and Elevators

CH82

Stripper 2.1cm wide,
elevator 1.2cm wide.
Overall length
5-1/2in (14cm).

CH80

Stripper 2.5cm wide,
elevator 1.6cm wide.
Overall length
8-3/4in (22.2cm).



CORYLLOS Raspatory

Double-ended.
Blades 0.8cm and 1.8cm.

CH90-1

Overall length
12in (30.5cm).

**OVERHOLT Periosteal Elevator**

Double-ended.

CH120

Blades
6.4mm and 7.1mm wide.
Overall length
11-1/2in (29.2cm).

**OVERHOLT Periosteal Elevator**

Double-ended.

CH130

Blades
6.4mm and 13mm wide.
Overall length
10-5/8in (27cm).

**CAMERON HAIGHT Periosteal Elevators**

Double-ended.
Round blade 12mm wide.
Straight edge 15mm wide.

CH145

Overall length
9-1/4in (23.5cm).

CH147

Overall length
13in (33cm).

**ALEXANDER-FARABEUF Costal Periosteotome**

Double-ended.

CH150

Blades 12mm x 21mm wide.
Overall length
8-1/4in (20.6cm).

**SEDILLOT Periosteal Elevator****CH165**

Broad rounded blade
19mm wide.
Overall length
7-1/2in (19.1cm).



Rib Shears

Costotome
Rippenscheren
Cesoie costali
Esquileo de costillas
Tesoura torácica
肋骨剪刀

POTTS Extra Small Rib Shears

CH290

Curved,
serrated blades.
Overall length
6in (15.2cm).



RUSKIN Bone Splitting Forceps

Double-action.
Overall length
7-3/8in (18.7cm).

◆ **CH470**

Straight blades.

CH471

Angled blades.



GLUCK Rib Shears

CH315

Overall length
7-3/4in (19.6cm).



GIERTZ-SHOEMAKER Rib Shears

◆ **CH330**

Overall length
9-3/4in (24.8cm).



SHOEMAKER Rib Shears

CH431

Overall length
9in (22.9cm).




Costotome
Rippenscheren
Cesoie costali
Esquileos de costillas
Tesoura torácica
肋骨剪刀

SAUERBRUCH Rib Shears

CH440

Overall length
10in (25.4cm).



CORYLLOS Rib Shears

Protected points.
Overall length
14-1/4in (36.2cm).

CH360

Right.

CH362

Left.



BETHUNE Rib Shears

CH340

Overall length
14in (35.6cm).



CORYLLOS-BETHUNE Rib Shears

Protected points.
Overall length
14-1/2in (36.8cm).

CH342

Right.

CH344

Left.



COOLEY Rib Shears

Angled cutting surface.

CH350

Heavy,
curved blades.
Overall length
14in (35.6cm).



TUDOR-EDWARDS Costotome

☐ **CH450**

Overall length
9-1/4in (23.5cm).




Rongeur Forceps

Pince-gouge
Rongeur-Zangen
Pinza ossivora
Fórceps de gúbias
Pinça saca-bocado
骨鉗子


SAUERBRUCH Rib Rongeur Forceps

Double-action. Square jaws.
Overall length 12in (30.5cm).




CH560

Straight jaws,
14mm wide.




CH564

Straight jaws,
20mm wide.



CH584

Curved jaws,
20mm wide.



DALE Rib Rongeur

Jaws 15mm x 3mm wide.



CH590

Overall length
14-1/2in (36.8cm).



LEBSCHKE Sternum Knife

CH680

Overall length
10in (25.4cm).



Sternal Perforating Awl

Aluminum Handle.
Needle length
1-3/8in (34mm).

CH691

Overall length
5-1/2 inch (14cm).



SAUERBRUCH Retractors

Flat blade is slightly curved on end.
Overall length 9-1/2in (24.1cm).

CH810

Blade
1-3/4in x 5/8in
(4.4cm x 1.6cm) wide.

CH805

Blade
2-1/4in x 7/8in
(5.7cm x 2.2cm) wide.

CH800

Blade
2-7/8in x 3/4in
(7.3cm x 1.9cm) wide.



SEMB Lung Retractors

Overall length 9-1/4in (24.8cm).



CH830

Blade
1-1/8in x 3/4in
(2.8cm x 1.9cm) wide.



CH832

Blade
1-5/8in x 1in
(4.1cm x 2.4cm) wide.



ALLISON Lung Retractors

CH840

Blade 2-1/8in (5.4cm) wide.
Overall length
12-1/4in (31.1cm).

CH841

Blade 1-1/2in (3.8cm) wide.
Overall length
10-3/8in (26.4cm).



DAVIDSON Scapular Retractors

CH905

Blade 3in x 2in (7.6cm x 5.1cm) wide.

CH900

Blade 3in x 3-1/2in (7.6cm x 8.9cm) wide.



Retractors

Écarteurs
Wundspreizer/Retraktoren
Divaricatori
Separadores
Retradores
リトラクター

PARSONNET Cardiovascular Retractors

Self-retaining,
light spring steel retractor
with 3 x 3 sharp prongs.

CH5751

Overall length
1-1/4in (3.2cm).

CH5752

Overall length
1-1/2in (3.8cm).

◆ CH5750

Overall length
2in (5.1cm).




COOLEY Atrial Retractors

Right. Serrated blade. Overall length 9-1/4in
(23.5cm).

CH8119

Blade 7/8in x 3/4in
(2.2cm x 1.9cm) wide.



CH8121

Blade 1-1/8in x 3/4in (2.9cm x 1.9cm) wide.

CH8123

Blade 1-1/2in x 1-1/2in (3.8cm x 3.8cm) wide.

CH8124

Blade 1-7/8in x 1-3/4in (4.8cm x 4.4cm) wide.

CH8120

Blade 2in x 1-3/4in (5.1cm x 4.4cm) wide.



ROCHESTER Atrial Septum Retractors

CH8152

Blade 5mm.
Overall length
5-1/2in (14cm).

CH8153

Blade 7.5mm.
Overall length
5-1/2in (14cm).

CH8154

Blade 10mm.
Overall length
5-7/8in (14.9cm).

CH8155

Blade 12mm.
Overall length
5-3/4in (14.6cm).

CH8156

Blade 15mm.
Overall length
6in (15.2cm).




COOLEY Atrial Retractors

Left. Serrated blade.

CH8127

Blade 3/4in x 3/4in
(1.9cm x 1.9cm) wide.
Overall length
9in (23cm).



CH8126

Blade 3/4in x 1-1/2in
(1.9cm x 3.8cm) wide.
Overall length
9-1/2in (24.1cm).

CH8125

Blade 3/4in x 1-3/4in
(1.9cm x 4.4cm) wide.
Overall length
9-1/2in (24.1cm).



CRAWFORD Aortic Retractors

Overall length
7-1/2in (19cm).

CH8110

Small.
Blade
15mm x 15mm.

CH8112

Large.
Blade
20mm x 22mm wide.


COOLEY Aorta Retractors

Double ended. Blades have
fine teeth along edges.
Blade 1/2in (1.3cm) wide.
Overall length
8in (20.3cm).

CH8130

Blade length
1-1/2in (3.8cm) and
5/8in (1.6cm).

CH8131

Blade length
1-3/4in (4.4cm) and
7/8in (2.2cm).


REVUELTA-GARCIA Atrial Retractors


Overall length
10in (25.4cm).

CH8114

Blade
1in (2.8cm) wide.



ROSS Aortic Valve Retractors

Overall length 10in (25.4cm)




CH8105-001

No. 1 blade, 1.5cm x 1.6cm wide.




CH8105-002

No. 2 blade, 1.3cm x 2.3cm wide.




CH8105-003

No. 3 blade, 1.8mm x 2.5mm wide.




CH8105-004

No. 4 blade, 1.3cm x 2.9cm wide.



CH8105-005

No. 5 blade, 2.9cm x 2.9cm wide.



CH8105-006

No. 6 blade, 3.5cm x 2.1cm wide.




Retractors/Forceps

Écarteurs/Forceps
Wundspreizer/Retraktoren/Zangen
Divicatori/Pinze
Separadores/Fórceps
Retradores/Fórceps
リトラクター/鉗子

ROCHESTER Rake Retractor

CH8158

3 sharp prongs.
8mm wide.
Overall length
6-3/4in (17.1cm).



CH8270

Holder 5in (12.7cm) wide. Overall length
5-1/4in (13.5cm).



CASTANEDA Suture Tag Forceps


CH6805

Smooth jaw.
Overall length
4in (10.2cm).




KARMODY Vascular Spring Retractors

Blades have serrated edges.




CH8669

Small. Overall length
3-1/8in (7.9cm).



CH8670



Large. Overall length
4-3/4in (12.1cm).



Pincers-gouges/Râteaux/Écarteurs
Rongeur/Spreizer/Wundspreizer/Retraktoren
Pinze ossivore/Divaricatori/Divaricatori
Gubias/Espátulas/Separadores
Saca-bocados/Afastadores/Retratores
骨鉗子/スプレッダー/リトラクター

BAILEY Rongeurs

Jaws 3mm wide.
Shaft length
4-15/16in (12.5cm).



CH1400

Straight.



CH1402

Jaws angled down 25°.

FINOCHIETTO Pediatric Rib Spreader

Curved blade arms.



◆ **CH1210-1**

Blades 1-1/4in x 1-3/4in (3.2cm x 4.4cm) wide.
Arm length 5-1/2in (14cm). Spread 6in (15.2cm).

FINOCHIETTO Pediatric Rib Spreader

Aluminum arms and blades.



CH1211

Blades 1-1/4in x 1-3/4in (3.2cm x 4.4cm) wide.
Arm length 4-5/16in (11cm).
Spread 6-3/8in (16.2cm).

COOLEY Rib Retractors

Pediatric. Aluminum arms and blades.



CH1249

Blades 16mm x 25mm wide.
Arm length 4in (10.2cm).
Spread 4-3/4in (12.1cm).

CH1248

Blades 35mm x 35mm wide.
Arm length 5-3/4in (14.6cm).
Spread 7in (17.8cm).

WILSON Rib Spreader

Aluminum arms and blades.



CH1250

Blades. 26mm x 30mm wide.
Arm length 4in (10.2cm).
Spread 3-7/8in (9.8cm).

Retractors/Spreaders

Écarteurs/Râteaux
Wundspreizer/Retraktoren/Spreizer
Divaricatori/Divaricatori
Separadores/Espátulas
Retratores/Afastadores
リトラクター/スプレッダー

KRASKY Retractor

Curved blade with six sharp prongs.
Blade 2in x 1in
(5.1cm x 2.5cm) wide.

CH860

Overall length
11in (27.9cm).



WEINBERG Rib Spreader

Aluminum arms and blades.



CH1262

Blades. 12mm x 19mm wide. Arm length
3-1/2in (8.9cm). Spread 2-7/8in (7.3cm).

HAIGHT Rib Retractor



CH1251

Large. Blades 1-1/4in x 1-1/8in (3.2cm x 2.7cm) wide. Arm length 3-1/2in (8.9cm).
Spread 4-1/2in (11.4cm).

FINOCHIETTO Rib Spreaders

Aluminum arms and blades.



CH1260

Blades. 10mm x 15mm wide. Arm length
2-3/4in (7cm). Spread 2-3/4in (7cm).

◆ CH1261

Blades. 6mm x 6mm wide. Arm length 1-3/8in
(3.5cm). Spread 1-3/8in (3.5cm).

BURFORD Rib Spreader

Aluminum arms and blades.



CH1243

Two blade sets: 1-1/4in x 1-1/8in (3.2cm x 2.9cm) wide and 1-3/4in x 1-1/8in (4.5cm x 2.9cm) wide. Arm length 6in (15.2cm). Spread
5-1/2in (14cm).

RIENHOFF-FINOCHIETTO Rib Spreader




CH1140-1

Blades 3/4in x 1in (1.9cm x 2.5cm) wide. Arm length 5-5/8in (14.3cm). Spread 6in (15.2cm).

LEMMON Sternal Spreader

Each arm has two blades.



CH1230

Blades 3.3cm x 3.1cm wide. Arm lengths 7-7/8in (20cm) and 6-5/8in (16.8cm). Spread 8in (20.3cm).

COOLEY Sternal Retractor

Extra small. Aluminum blades and arms.



CH6950-007

Blades 3/8in x 5/8in (10mm x 16mm) wide. Arm length 2-1/8in (5.4cm). Spread 3-1/8in (7.9cm).

COOLEY Sternotomy Retractor

Small. Aluminum arms and blades. With stabilizing pins in blades.




CH1289

Blades 14mm x 19mm wide. Arm length 3-3/8in (8.6cm). Spread 4in (10.2cm).

COOLEY Sternotomy Retractors

Aluminum arms and blades.



CH1287

Blades 20mm x 29mm wide. Arm length 3-1/8in (7.9cm). Spread 4-5/8in (11.7cm).

CH1288

Blades 14mm x 16mm wide. Arm length 2-3/4in (7cm). Spread 3-3/8in (8.6cm).

Retractors/Spreaders

Écarteurs/Râteaux
Wundspreizer/Retraktoren/Spreizer
Divaricatori/Divaricatori
Separadores/Espátulas
Retratores/Afastadores
リトラクター/スプレッダー

HARRINGTON Retractors

Overall length 14-1/2in (36.8cm).

CH1008

Large blade 7-1/2in x 4-1/2in
(19.1cm x 11.4cm) wide.

CH1009

Extra-large blade 7-1/2in x 6-1/2in
(19.1cm x 16.5cm) wide.



HARRINGTON-PEMBERTON Retractor

Blade 7-1/2in x 2-1/2in
(19.1cm x 6.4cm) wide,
thickened distal end.

CH1010

Overall length
14-1/2in (36.8cm).



HARRINGTON Splanchnic Retractors

Flat, deeply curved, springy blade with
thickened distal end.

Overall length
12-1/2in (31.7cm).

CH1000

Blade 5in x 1-1/2in
(12.7cm x 3.8cm) wide.


CH1005

Blade 5in x 2-1/2in
(12.7cm x 6.4cm) wide.



Titanium ANKENY Sternal Retractor

Includes six detachable swivel blades and one pair of springs on each arm.



CH1290

Blades.
1-1/4in x 1-3/8in wide
(3.2cm x 3.5cm).
Arm length 11-1/2in (29.2cm).
Spring length 3-1/4in (8.3cm).
Spread 7in (17.8cm).

CH1290-1

Blades only.
2in x 1-3/8in wide
(5.2cm x 3.5cm).
Six per set.

titanium

COOLEY-MERZ Sternum Retractor



CH1275-1

Extra-large.
Blades 1-3/8in x 4-3/8in
(3.5cm x 11cm) wide.
Arm length 8-1/2in (21.6cm).
Spread 9-3/8in (23.8cm).

COOLEY Sternotomy Retractors



CH1285-1

Blades 7/8in x 1-7/8in
(2.1cm x 4.6cm) wide.
Arm length 4-3/8in (11cm).
Spread 8-1/4in (21cm).

◆ CH1280-1

Blades.
1-3/8in x 2in (3.5cm x 5cm) wide.
Arm length 5-3/4in (14.5cm).
Spread 9-1/2in (24.1cm).

TUFFIER Rib Spreader




◆ **CH1100**

Blades 1-1/2in x 1-5/8in
(3.8cm x 4.1cm) wide.
Spread 7in (17.8cm).

TUFFIER Rib Retractor

Total opening 6-1/2in (16.5cm).



CH1101

Serrated blades 2in x 1-3/4in wide
(5.1cm x 4.4cm).
Blades 2in deep (5.1cm).

Spreaders

Râteaux
Spreizer
Divaricatori
Espátulas
Afastadores
スプレッダー

FINOCHIETTO Rib Spreader

Small. Straight arms.



CH1265

Blades 1/2in x 5/8in wide
(1.3cm x 1.6cm).
Arm length 2-1/16in (5.2cm).
Spread 2-2/8in (5.4cm).

FINOCHIETTO Rib Spreaders

Curved blade arms.
Arm length 7-1/4in (18.4cm).
Spread 8-1/2in (21.6cm).



◆ **CH1215-1**

Adult.
Blades 1-3/4in x 2-1/2in
(4.4cm x 6.4cm) wide.

CH1220-1

Large.
Blades 2-3/4in x 2-1/2in
(7cm x 6.4cm) wide.

FINOCHIETTO Rib Spreaders

Aluminum arms and blades.




CH1216

Adult.
Blades 1-1/2in x 2-3/8in (3.8cm x 6cm) wide.
Arm length 7-1/4in (18.4cm).
Spread 7-5/8in (19.4cm).

CH1221

Large.
Blades 3in x 2-3/8in (7.6cm x 6cm) wide.
Arm length 7in (17.8cm).
Spread 7-1/2in (19.1cm).

BURFORD Rib Spreader




◆ **CH1240-1**

Two small blades
1-3/4in x 2-1/2in
(4.4cm x 6.4cm) wide and
two large blades 3in x 2-1/2in
(7.6cm x 6.4cm) wide.
Spread 7-1/2in (19.1cm).

BURFORD Rib Spreader

Aluminum arms and blades.



CH1242-1

Adult.
Two blade set: 2in x 2-1/4in
(5.1cm x 5.7cm) wide and
2-1/2in x 2-1/2in (6.4cm x 6.4cm) wide.
Arm length 9-5/16in (23.7cm).
Spread 7-7/8in (20cm).

BAILEY-GIBBON Rib Contractor

Large thumbscrew automatic ratchet holds blades in place.

CH1355

Long arms,
length 2in (5.1cm).
Overall length
6-1/2in (16.5cm).



LEMMON Sternal Approximator

Used in pairs to maintain close approximation of sternal edges.

CH1340

Prongs
1-1/16in (2.7cm).
Arm length 1-3/8in (3.5cm).
Overall length
7-1/8in (18.1cm).



BAILEY Rib Contractor

Large thumbscrew automatic ratchet holds blades in place.

◆ CH1350

Short arms,
length 1-3/8in (3.5cm).
Overall length
6-1/2in (16.5cm).



Customer support

V. Mueller products are sold only by authorized V. Mueller representatives. All instruments and equipment of our manufacture bear the V. Mueller trademark.

For more information, please contact your sales representative, send an email to **GMB-VM-Europe@carefusion.com** or visit **carefusion.com/vmueller**.

Due to changing surgical trends, some products described in this catalog may be discontinued or not available in your country.

Your V. Mueller sales representative may be able to suggest a suitable substitute if one is available.


Retractors

Écarteurs
Wundspreizer/Reraktoren
Divaricatori
Separadores
Retratores
リトラクター

HENLY Retractor Set

◆ CH1030

In addition to the retractor frame, the set includes three pair of side blades, each having four blunt prongs: 3in (7.6cm), 2in (5.1cm) and 1in (2.5cm). Also, four center blades: 3in (7.6cm), 1-1/4in (3.2cm), 1in (2.5cm) and 3/4in (1.9cm).



Cat. No.	Description	Set Quantity
CH1030-001	Center blade. 3/4in x 5/8in (1.9cm x 1.6cm) wide.	1
CH1030-002	Center blade. 1in x 5/8in (2.5cm x 1.6cm) wide.	1
CH1030-003	Center blade. 1-1/4in x 5/8in (3.2cm x 1.6cm) wide.	1
CH1030-007	Center blade. 3in x 5/8in (7.6cm x 1.6cm) wide.	1
CH1030-004	Side blades. 1in x 1in (2.5cm x 2.5cm) wide.	Pair
CH1030-005	Side blades. 2in x 1in (5.1cm x 2.5cm) wide.	Pair
CH1030-006	Side blades. 3in x 7/8in (7.6cm x 2.2cm) wide.	Pair
CH1031	Frame only. Overall length 7in (17.8cm)	1

HENLY Retractors

Original-type. Modified version of the WEIT-LANER retractor allows for a third retractor blade to be attached.

◆ CH1130

Complete set includes HENLY retractor, overall length 6-1/2in (16.5cm), and three center blades (CH1131, CH1132, CH1133).

CH1131


Small.
3/4in x 5/8in wide (1.9cm x 1.6cm).

CH1132


Medium.
1in x 5/8in wide (2.5cm x 1.6cm).

CH1133

Large.
1-1/4in x 5/8in wide. (3.2cm x 1.6cm).



Dissecteurs/Crochets
Dissektoren/Haken
Dissettori/Uncini
Disseccionadores/Ganchos
Escalpelos/Ganchos
切開器/フック

Double-Ended Dissector**CH2605**

Overall length
7in (17.8cm).

OLIVEKRONA Double-Ended Dissectors**CH2615**

Light. Overall length
7-1/4in (18.4cm).

CH2616

Light. Overall length
9-5/8in (24.4cm).

CH2617

Heavy. Overall length
7-3/4in (19.7cm).

CH2618

Heavy. Overall length
9-1/2in (24.1cm).

CRILE Hook and Dissector

Matte and Mirror finish.

CH2600

Overall length
12-1/2in (31.8cm).

Villavicencia Vein Hooks

Overall length
5-3/4in (14.6cm).

CH2611

2mm hook.

CH2612

3mm hook.

CRILE Single Hook

Mirror and Satin finish.

CH2610

Overall length
12-1/2in (31.8cm).

Vascular Hooks

Overall length
9-1/2in (24.1cm).

CH8676

Angled 90°.

Scissors

Ciseaux
Scheren
Forbici
Tijeras
Tesoura
鉋

CASTROVIEJO Micro Scissors


Round handle.
Sharp fine tips.
Overall length
4-1/2in (11.4cm).

CH07420

Angled 25°.

CH07422

Angled 45°.



CH07424


Angled 60°.

CH07426

Angled 90°.

CH07428

Angled 125°.



Artery Scissors


Ball tip on lower blades.
Overall length
6-3/4in (17.2cm).

CH5661-001

Angled 25°.

CH5661-002

Angled 45°.



CH5661-004

Angled 90°.

CH5661-005

Angled 120°.



Artery Scissors

Overall length
6-3/4in (17.2cm).

CH5663-001

Blades angled 25°.

CH5663-002

Blades angled 45°.

CH5663-004

Blades angled 90°.




CH5663-005

Blades angled 120°.



Artery Scissors

Delicate.
Overall length
8in (20.3cm).



CH5668-002

Blades angled 45°.

CH5668-003

Blades angled 60°.




JACOBSON Angled Spring-Handled Scissors

Flat handle. Fine blades.
Overall length 6-1/2in (16.5cm).

CH07450

Angled 25° sharp.



CH07452

Angled 45° sharp.

CH07454

Angled 60° sharp.

CH07456

Angled 90° sharp.



CH07458

Angled 125° sharp.



JACOBSON Spring-Handled Scissors

Curved blades,
length 7/8in
(2.2cm) from mid-screw.
Mirror finish.

CH8545

Overall length
5-1/4in (13.3cm).

CH8550

Overall length
6-1/2in (16.5cm).



JACOBSON Micro Scissors

Flat handle.
Overall length
7in (17.8cm).

CH07470

Angled on flat blades 25°.

CH07472

Angled on flat blades 45°.

CH07274

Angled on flat blades 60°.

CH07276

Angled on flat blades 90°.

**JACOBSON Angled Spring-Handled Scissors**

Round handle. Fine blades.
Overall length 7in (17.8cm).

**CH07460**

Angled 25° sharp.

◆ **CH07462**


Angled 45° sharp.

CH07464

Angled 60° sharp.

CH07466

Angled 90° sharp.

◆ **CH07468**


Angled 125° sharp.

**JACOBSON Micro Scissors**


Sharp fine tips.
Overall length 7in (17.8cm).

CH07460-1

Blades angled 25°.

**CH07462-1**

Blades angled 45°.

**CH07464-1**

Blades angled 60°.

CH07466-1

Blades angled 90°.

CH07468-1

Blades angled 125°.

**JACOBSON Bayonet-Shaped Scissors**

Spring handles. Very delicate blades, angled 90° on-the-flat. Blade length 12mm from inside of curve to tips. Mirror and Satin finish. Overall length 7-3/8in (18.7cm).

**CH8555**


Sharp tips.

**CH8560**

Blunt tips.

**MILLS Circumflex Scissors**

Blades angled 125°.

**CH5706**

1mm ball on tip.
Overall length
7-1/4in (18.4cm).



Scissors

Ciseaux
Scheren
Forbici
Tijeras
Tesoura
鉋

YASARGIL Scissors

Bayonet-shaped.
Reverse angle cutting blades,
length 12mm.
Angled 125°.



CH5708

Overall length
7-1/4in (18.4cm).



LINDLEY Scissors

Delicate. Spring handled.
Blades angled 90°.

CH8360

Overall length
8in (20.3cm).



Customer support

V. Mueller products are sold only by authorized V. Mueller representatives. All instruments and equipment of our manufacture bear the V. Mueller trademark.


For more information, please contact your sales representative, send an email to **GMB-VM-Europe@carefusion.com** or visit **carefusion.com/vmueller**.

Due to changing surgical trends, some products described in this catalog may be discontinued or not available in your country.

Your V. Mueller sales representative may be able to suggest a suitable substitute if one is available.


DIETHRICH Artery Scissors

V. Mueller offers the full family of DIETHRICH scissors that come in five different blade angles, available in 25°, 45°, 60°, 90°, and 125° blade angles. Scissors have tapered shafts with thin, narrow blades provide delicate, lightweight feel; finer than the POTTS-SMITH style. Satin finish prevents glare. Blade length 3/4in (19mm). Overall length 7-1/4in (18.4cm).




CH5725

Blades angled 45°.



CH5730

Blades angled 60°.




CH5735

Blades angled 90°.



CH5720

Blades angled 25°.



CH5740

Blades angled 125°.
With stabilizer bar.



DIETHRICH Artery Scissors

Extra delicate blades.
Overall length
7in (17.8cm).

CH5720-001


Blades angled 25°.

CH5725-001

Blades angled 45°.

CH5730-001

Blades angled 60°.

**CH5735-001**

Blades angled 90°.

CH5740-001

Blades angled 120°.

**Super Cut DIETHRICH Artery Scissors**

Sharp/Sharp. One serrated blade. Mirror finish.

CH5666-03S


Angled 90°.
Overall length 6-3/4in (17.1cm).

CH5666-02S

Angled 60°.
Overall length 7in (17.8cm).

CH5666-01S

Angled 25°.
Overall length 7-1/4in (18.4cm).

**CH5666-S**

Angled 45°.
Overall length 7in (17.8cm).

CH5705-S


With stabilizing bar.
Angled 125°. Overall length 7in (17.8cm).

**POTTS-SMITH Scissors**

Blade length 5/8in (1.5cm).
Overall length 5in (12.7cm).

CH5651-001

Blades angled 25°.

**CH5651-002**


Blades angled 60°.

CH5651-003

Blades angled 90°.

**POTTS-SMITH Scissors****CH5651**

Short. Angled 45°.
Blade length 1/2in (1.3cm).
Overall length
5-1/4in (13.3cm).



Scissors

Ciseaux
Scheren
Forbici
Tijeras
Tesoura
鉋

POTTS-Type Cardiovascular Scissors

Tapered shanks. Thin, narrow, micro blunt blades, 10mm.

CH5671-001

Angled 25°. Overall length 7in (17.8cm).

CH5671-002

Angled 45°. Overall length 7in (17.8cm).

CH5671-003

Angled 60°. Overall length 6-7/8in (17.5cm).

CH5671-004

Angled 90°. Overall length 6-3/4in (17.1cm).

CH5671-005

Angled 120°. Overall length 6-3/4in (17.1cm).



POTTS-SMITH Scissors

Regular blades, length 7/8in (22mm). Mirror finish. Overall length 7-1/4in (18.4cm).

CH5665

Angled 25°.

◆ CH5662

Angled 45°.

CH5660

Angled 60°.



POTTS-SMITH Scissors

Delicate. Blade length 1/2in (13mm).

CH5666-001

Angled 25°. Overall length 7-1/4in (18.2cm).

◆ CH5666


Angled 45°. Overall length 7in (17.8cm).

CH5666-002

Angled 60°. Overall length 7in (17.8cm).

CH5666-003

Angled 90°. Overall length 6 3/4in (17.1cm).




POTTS-SMITH Scissors

Mirror finish.

CH5652

Short. Blades angled 45°, length 3/4in (19cm). Overall length 5-3/4in (14.6cm).



POTTS-SMITH Scissors

Blades angled-on-flat 50°. Mirror finish.

CH5670

Blade length 1in (25mm). Overall length 7-3/4in (19.7cm).



Super Cut POTTS-SMITH Scissors

Mirror finish.
Overall length
7in (17.8cm).

CH5665-S

Angled 25°.

CH5662-S


Angled 45°.

CH5660-S

Angled 60°.

**POTTS-SMITH Reverse Scissors**

Blades angled 130°,
length 1/2in (1.3cm).
Overall length
7in (17.8cm).



◆ **CH5704**

Without stabilizing bar.

CH5705

With stabilizing bar.

**COOLEY Reverse Cut Scissors**

Mirror finish.

CH5700

Blades angled 120°.
Overall length
6-3/4in (17.1cm).

**DeBAKEY-POTTS Scissors**

Blade length 1in (25mm).
Mirror finish.
Blunt/blunt.
Overall length
6-1/2in (16.5cm).

CH5667-001

Blades angled 25°.

CH5667-002

Blades angled 45°.

CH5667-003

Blades angled 60°.

**DeBAKEY Vascular Scissors**

Blade length 1in (25mm).
Mirror finish.

CH5669-001

Blades angled 25°.
Overall length
9in (22.9cm).

CH5669-002


Blades angled 45°.
Overall length
9in (22.9cm).

CH5669-003

Blades angled 60°.
Overall length
9in (22.9cm).

CH5670-001

Blades angled 25°.
Overall length
11-1/4in (28.6cm).

**COOLEY Vascular Scissors**

Overall length
7-1/4in (18.4cm).
Mirror finish.
Blades angled 45°.

CH5656

1mm probe
point on tip.

CH5657

1.5mm probe
point on tip.




Scissors

Ciseaux
Scheren
Forbici
Tijeras
Tesoura
鋏

COOLEY Arteriotomy Scissors

Blade angle 30°,
length 1/2in (12mm).
Mirror finish.



CH5672

Overall length
7-1/4in (18.4cm).

CH5673

Overall length
9in (22.9cm).

CH5674

Overall length
10-3/4in (27.3cm).




MORSE Scissors

Sharply pointed blades are
externally sharpened.
Bottom blade serrated.
Mirror finish.

CH5695

Curved.
Overall length
7in (17.8cm).



FAVALORO Curved Shank Coronary Scissors

Blades angled 130°.

CH5653

Blade length
7/16in (11mm).
Overall length
6-3/4in (17.1cm).



POTTS Tenotomy Scissors

Overall length
5-1/2in (14cm).

◆ **CH5674-001**

Curved.

CH5674-002

Straight.



POTTS Tenotomy Scissors

Curved.
Mirror finish.

◆ **CH5675**

Overall length
6in (15.2cm).

◆ **CH5676**

Overall length
7in (17.8cm).

CH5677

Overall length
8in (20.3cm).



LILEHEI-POTTS Tenotomy Scissors

Overall length
7in (17.8cm).

CH0758

Straight.

CH0764

Curved.



Super Cut JAMISON Scissors

Curved.
Mirror finish.

CH2035-S

Overall length
5-1/2in (14cm).

CH2036-S


Overall length
7in (17.8cm).

**JAMISON-METZENBAUM Scissors**

Curved.
Tenotomy type tips.

CH5677-9

Overall length
9in (22.8cm).

**CHURCH Scissors**

Light, long scissors.
Beveled, MAYO-type
dissecting blades.
Blunt.

CH2240

Straight.
Overall length
6-3/4in (17.1cm).

CH2260


Straight.
Overall length
10in (25.4cm).

CH2250

Curved.
Overall length
6-3/4in (17.1cm).

CH2255

Curved.
Overall length
10in (25.4cm).

**POTTS-SMITH Dissecting Scissors**

Curved. Mirror finish.

CH5681

Overall length
5-1/2in (14cm).

◆ **CH5680**

Overall length
7-1/2in (19.1cm).

CH5679

Overall length
9in (22.9cm).

CH5678

Overall length
11-1/4in (28.6cm).

**NELSON Lobectomy Scissors**

Rounded, curved blades.

CH2026

Overall length
10-1/4in (26cm).

CH2027

Overall length
14-1/4in (36.2cm).

**COOLEY Cardiovascular
"My" Scissors**

MAYO-type tips.
Curved.
Mirror finish.

◆ **CH2052**

Overall length
7-1/2in (19.1cm).

◆ **CH2053**

Overall length
9in (22.9cm).



Scissors

Ciseaux
Scheren
Forbici
Tijeras
Tesoura
鉋

SWEET Esophageal Scissors

Blades angled on flat blunt tips.

CH2060

Overall length
8-1/2in (21.6cm).



JORGENSEN Scissors

A modification of the MAYO dissecting scissors with an exaggerated curve. Mirror finish.

◆ **CH2070**

Overall length
8-1/2in (21.6cm).



LITWIN Angled Scissors

Blades angled 45°. Mirror finish.

CH5655

Overall length
7-1/4in (18.4cm).



LITWAK General Utility Scissors

Curved. Blunt/Blunt. Mirror finish.

CH5688

Overall length
4-5/8in (11.7cm).



CRAFOORD Lobectomy Scissors

Curved heavier blade. Blunt/Blunt.

CH2040

Overall length
12in (30.5cm).
Mirror and Satin finish.



Cardiovascular Scissors

Very fine and delicate. Small ball tip on one blade. Straight. Mirror finish.

CH5690

Overall length
4-3/4in (12.1cm).



DeBAKEY Endarterectomy Scissors

Mirror finish.



CH5697

Curved. Blunt/Blunt.
Overall length
7in (17.8cm).



Super Cut Endarterectomy Scissors

Curved.
Blunt/Blunt.
One serrated blade.
Mirror finish.

CH5698-S

Overall length
7in (17.8cm).



TOENNIS Dissecting Scissors

Delicate. Mirror finish.

CH5682

Straight.
Overall length
7-1/4in (18.4cm).

CH5684

Curved.
Overall length
7-1/8 (18.1cm).



DeMARTEL Vascular Scissors

Slim and delicate.
Blunt/Blunt.
Blade length
from mid-screw
1-1/8in (2.9cm).
Mirror finish.

CH5693


Straight.
Overall length
7-3/4in (19.7cm).

CH5692

Curved.
Overall length
7-3/4in (19.7cm).

CH5692-001

Curved.
Overall length
6-3/4in (17.1cm).



THOMPSON Valve Scissors

Fully curved blades,
one serrated blade.

CH0812

Overall length
8-1/2in (21.6cm).

CH0814

Overall length
9-7/8in (25cm).



LITWAK Cardiovascular Scissors

Angled-on-side.
Overall length
11in (27.9cm).

CH5686

Slightly curved blades.

CH5687

Fully curved blades.




Scissors/Forceps

Ciseaux/Forceps
Scheren/Zangen
Forbici/Pinze
Tijeras/Fórceps
Tesoura/Fórceps
鉗/鉗子

Valve Leaflet Excision Scissors

Spring handle.
Blades angled
45° on-the-flat.




CH8689

Overall length
5-3/4in (14.6cm).



Arteriotomy (Finger POTTS) Scissors

Spring handle.
Blades angled 60°.



◆ **CH8690**

Overall length
5-3/8in (13.7cm).



COOLEY Scissors

Delicate, curved METZENBAUM-type blades
Gold shanks for identification.

CH6950-008

Overall length
5-3/4in (14.6cm).



Customer support

V. Mueller products are sold only by authorized V. Mueller representatives. All instruments and equipment of our manufacture bear the V. Mueller trademark.

For more information, please contact your sales representative, send an email to **GMB-VM-Europe@carefusion.com** or visit **carefusion.com/vmueller**.

Due to changing surgical trends, some products described in this catalog may be discontinued or not available in your country.


Your V. Mueller sales representative may be able to suggest a suitable substitute if one is available.

LEATHER KARMODY In Situ Valve Scissors

Blunt-tipped, serrated blades. Alligator-style.

**CH8675**Working length
6-3/4in (17.2cm)**MAYO-HARRINGTON
Dissecting Scissors**

Mirror finish.

CH2085Straight.
Overall length
9in (22.9cm).**CH2080**Straight.
Overall length
11in (27.9cm).◆ **CH2095**Curved.
Overall length
9in (22.9cm).**CH2090**Curved.
Overall length
11in (27.9cm).**METZENBAUM Dissecting
Scissors**

Blunt/Blunt. Mirror and Satin finish.

CH2034-001Straight. Overall length
8in (20.3cm).**CH2034**Curved. Overall length
8in (20.3cm).◆ **CH2032**Curved. Overall length
9in (22.9cm).◆ **CH2030**Curved. Overall length
11-3/16in (28.4cm).**NELSON Lung
Dissecting Scissors**Blunt/Blunt.
Mirror and Satin finish.
Overall length
9in (22.9cm).**CH2010**

Curved.

CH2015

Straight.



Scissors

Ciseaux
Scheren
Forbici
Tijeras
Tesoura
鉋/鉋子

Super Cut METZENBAUM Dissecting Scissors

*Curved. One serrated blade.
Mirror and Satin finish.*

CH2032-S

*Overall length
9in (22.9cm).*

CH2025-S

*Overall length
11-3/4in (29.8cm).*



Super Cut MAYO Dissecting Scissors

Mirror finish.

CH2086-S

*Straight. Overall length
9in (22.9cm).*

CH2096-S

*Curved.
Overall length
9in (22.9cm).*



Super Cut MAYO-HARRINGTON Scissors

*Curved. One serrated blade.
Mirror and Satin finish.*

CH2075-S

*Overall length
12in (30.5cm).*



Vital® METZENBAUM Dissecting Scissors

*Delicate. Sharp/Sharp.
Stellite® cutting edges.*

CH2007-002

*Straight.
Overall length
9-1/4in (23.5cm).*

CH2006-002

*Curved.
Overall length
9-1/4in (23.5cm).*



Vital® NELSON Dissecting Scissors

Curved. Blunt/Blunt. Stellite® cutting edges.
Overall length 11in (27.9cm).

CH2020*Regular.***CH2020-001***Delicate.***Vital® MAYO Dissecting Scissors**

Beveled blades. Stellite® cutting edges.

CH2086

Straight.
Overall length
9in (22.9cm).

◆ **CH2096**

Curved.
Overall length
9in (22.9cm).

CH2088

Straight.
Overall length
14in (35.6cm).

CH2098

Curved.
Overall length
14in (35.6cm).

**Vital® METZENBAUM Dissecting Scissors**

Blunt/Blunt. Stellite® cutting edges.

CH2007

Straight.
Overall length
9-1/4in (23.5cm).

◆ **CH2006**


Curved.
Overall length
9-1/4in (23.5cm).

CH2009

Curved.
Overall length
14-1/2in (36.8cm).

CH2008

Curved.
Overall length
12in (30.5cm).

**Vital® METZENBAUM Dissecting Scissors**

Delicate. Blunt/Blunt. Stellite® cutting edges.

CH2007-001

Straight.
Overall length
9in (22.9cm).

◆ **CH2006-001**

Curved.
Overall length
9in (22.9cm).



Scissors/Wire Twisters/Needle Holders

Ciseaux/Cercleurs/Porte-aiguilles
Scheren/Drahtdreher/Nadelhalter
Forbici/Torcifilo/Porta-ago
Tijeras/Herramientas para enrollar alambres/Sostenes de agujas
Tesoura/Enroladores de fio/Porta-agulhas
鉋/ワイヤツイスター/針ホルダー

DIETHRICH Cardiovascular Scissors

Serrated blades.
Overall length
9in (22.9cm).

CH5741

Slightly curved blades.

CH5742

Fully curved blades.



SATINSKY Scissors

Fully curved blades,
slightly curved handle.
Mirror finish.

CH5650

Blade length
1-7/8in (4.8cm).
Overall length
9-3/4in (24.8cm).



Vital® RUBIO Wire Twister

Tungsten Carbide inserts.
3,600 teeth per square inch.
Overall length
5-1/4in (13.3cm).



◆ **CH2535**



Vital® Wire Twisters

Rounded jaws. Tungsten Carbide inserts.
3,600 teeth per square inch.
Overall length
6in (15.2cm).




◆ **CH2525**

Delicate.



CH2530

Heavy.



Vital® COOLEY-BAUMGARTEN Wire Twister

Heavy, blunt jaws. Tungsten Carbide inserts.
3,600 teeth per square inch.



◆ **CH2520**

Overall length
8in (20.3cm).



WANGENSTEEN Needle Holder

Straight jaw.



CH2432

Overall length
10-3/4in (27.3cm).




STILLE-FRENCH Cardiovascular Needle Holder

Delicate.



CH2398

Overall length
6-1/2in (16.5cm).



CRILE-MURRAY Needle Holder

Straight, cross-serrated jaws.



CH2400

Overall length
7-1/4in (18.4cm).



Vital® ROCHESTER Needle Holder

Light and delicate. Straight jaws.
Tungsten Carbide inserts.
3,600 teeth
per square inch.



◆ **CH2417**

Overall length
7-1/8in (18.1cm).



Vital® COOLEY Microvascular Needle Holders

Very delicate, straight jaws. Indented shanks.
Tungsten Carbide inserts.



CH2406

Smooth inserts.
Overall length
6in (15.2cm).



CH2408

Smooth inserts.
Overall length
8in (20.3cm).



CH2407-001

Serrated inserts.
3,600 teeth
per square inch.
Overall length
6-3/4in (17.1cm).



CH2407

Smooth inserts.
Overall length
6-3/4in (17.1cm).



CH2406-001

Serrated inserts.
3,600 teeth
per square inch.
Overall length
6in (15.2cm).



CH2408-001

Serrated inserts.
3,600 teeth
per square inch.
Overall length
8in (20.3cm).




Needle Holders

Porte-aiguilles
Nadelhalter
Porta-ago
Sostenes de agujas
Porta-agulhas
針ホルダー

Vital® COOLEY Intracardiac Needle Holders

Straight, narrow jaw. Tungsten Carbide inserts.
3,600 teeth per square inch.



CH2492

Overall length
7in (17.8cm).

CH2496

Overall length
9in (22.9cm).

CH2494


Overall length
8in (20.3cm).

CH2498

Overall length
10-3/8in (26.3cm).

Vital® Microvascular Needle Holders

Straight jaws, serrated. Tungsten Carbide inserts.
3,600 teeth per square inch.



CH2443

Overall length
7-1/8in (18.1cm).

CH2446

Overall length
10-1/4in (26cm).

CH2444


Overall length
8in (20.3cm).

CH2445

Overall length
9in (22.9cm).

Vital® RYDER Needle Holders

Narrow jaws. Tungsten Carbide inserts. Gold Plate and Satin finish.



CH2500

Smooth inserts.
Overall length
5in (12.7cm).

CH2501

Smooth inserts.
Overall length
6in (15.2cm).

CH2502

Smooth inserts.
Overall length
7in (17.8cm).

CH2503

Smooth inserts.
Overall length
8in (20.3cm).

CH2504

Smooth inserts.
Overall length
9in (22.9cm).

◆ **CH2506**

Serrated inserts.
3,600 teeth per
square inch.
Overall length
5in (12.7cm).

◆ **CH2507**

Serrated inserts.
3,600 teeth per
square inch.
Overall length
6in (15.2cm).

◆ **CH2508**

Serrated inserts.
3,600 teeth per
square inch.
Overall length
7in (17.8cm).

CH2509

Serrated inserts.
3,600 teeth per
square inch.
Overall length
8in (20.3cm).

CH2510

Serrated inserts.
3,600 teeth per
square inch.
Overall length
9in (22.9cm).

Vital® DeBAKEY Cardiovascular Needle Holders

Delicate patterns. Straight.
 Tungsten Carbide inserts.
 3,600 teeth per square inch.

◆ **CH2481**

Overall length
 7-1/8in (18.1cm).

**CH2483**


Overall length
 10-3/8in (26.3cm).

◆ **CH2482**

Overall length
 9in (22.9cm).

◆ **CH2484**

Overall length
 11-3/4in (29.8cm).

**Vital® Microsurgery Needle Holders**

Spring handle. Smooth Tungsten Carbide inserts.

**CH2448**

Curved jaws. Flat lock.
 Overall length
 5-5/8in (14.3cm).

**CH2452**

Straight jaws. Flat lock.
 Overall length
 5-3/4in (14.6cm).

CH2454

Straight jaws. Without lock.
 Overall length
 5-3/4in (14.6cm).

◆ **CH2457**

Straight jaws. Flat lock.
 Overall length
 7-1/4in (18.4cm).

**Vital® SPENCER Cardiovascular Needle Holder**

Spring handle.
 Flat lock. Smooth
 Tungsten Carbide inserts.

**CH8800**

Overall length
 7in (17.8cm).

**Vital® CASTROVIEJO Needle Holders**

Spring handle.
 Flat lock. Straight jaws.
 Tungsten Carbide inserts.
 3,600 teeth per square inch.

◆ **CH2456**

Overall length
 7-1/4in (18.4cm).

◆ **CH8589**

Overall length
 8-3/4in (22.2cm).

**Vital® JACOBSON Needle Holder**

Straight jaw.
 Smooth Tungsten
 Carbide inserts.

**CH8591**

Round lock.
 Overall length
 8-3/4in (22.2cm).



Needle Holders

Porte-aiguilles
Nadelhalter
Porta-ago
Sostenes de agujas
Porta-agulhas
針ホルダー

JACOBSON Needle Holders

Spring handle.
Round lock.
Straight, smooth jaws.



CH8592

Overall length
7in (17.8cm).



CH8590

Overall length
8-3/4in (22.2cm).



Vital® JACOBSON Spring-Handled Needle Holder

Straight jaws.
Smooth Tungsten
Carbide inserts.



CH8591-001

Flat lock.



JACOBSON Spring-Handled Needle Holders



CH8590-001

Overall length
8-5/8in (22cm).



CH8592-001

Overall length
7in (17.8cm).



Vital® JACOBSON Heavy-Style Needle Holders

Round handle for large needles.
Tungsten Carbide inserts.
3,600 teeth per square inch.
With lock.

CH06310

Overall length
7-1/2in (19cm).



◆ **CH06330**

Overall length
9in (22.9cm).

◆ **CH06340**

Overall length
10in (25.4cm).



Delicate Touch® CASTROVIEJO Micro Needle Holders

Round handle. Overall length
5-3/4in (14.6cm).



CH0650

Straight.
With lock.

CH0652

Straight.
Without lock.



CH0654

Curved.
With lock.

CH0656

Curved.
Without lock.



Delicate Touch® CASTROVIEJO Micro Needle Holders

Flat handle. Overall length
5-3/4in (14.6cm).



CH0658

Straight.
With lock.

CH0660

Straight.
Without lock.



CH0662

Curved.
With lock.

CH0664

Curved.
Without lock.



Delicate Touch® JACOBSON Micro Needle Holders

Flat handle.

CH0670

 Straight. With lock.
 Overall length
 6-3/4in (17.1cm).

CH0672

 Straight. Without lock.
 Overall length
 6-3/4in (17.1cm).

CH0674

 Curved. With lock.
 Overall length
 6-3/4in (17.1cm).

CH0676

 Curved. Without lock.
 Overall length
 6-3/4in (17.1cm).

CH0690

 Straight. With lock.
 Overall length
 7-1/4in (18.4cm).

CH0692

 Straight. Without lock.
 Overall length
 7-1/4in (18.4cm).

CH0694

 Curved. With lock.
 Overall length
 7-1/4in (18.4cm).

CH0696

 Curved. Without lock.
 Overall length
 7-1/4in (18.4cm).

Delicate Touch® JACOBSON Micro Needle Holders

Round handle. Mirror and Satin Finish.

 ◆ **CH0682**

 Straight. With lock.
 Overall length
 7-1/4in (18.4cm).

CH0684

 Straight. Without lock.
 Overall length
 7-1/4in (18.4cm).

 ◆ **CH06121**

 Straight. With lock.
 Overall length
 8-1/4in (21cm).

CH06123

 Straight. Without lock.
 Overall length
 8-1/4in (21cm).

CH0686

 Curved. With lock.
 Overall length
 7-1/4in (18.4cm).

CH0688

 Curved without lock.
 Overall length
 7-1/4in (18.4cm).

CH06125

 Curved. With lock.
 Overall length
 8-1/4in (21cm).

CH06127

 Curved. Without lock.
 Overall length
 8-1/4in (21cm).


Needle Holders

Porte-aiguilles
Nadelhalter
Porta-ago
Sostenes de agujas
Porta-agulhas
針ホルダー

Delicate Touch® JACOBSON Micro Needle Holders

Dust jaw surface for use with 5-0 suture and smaller. Straight. Overall length 8in (20.3cm). Flat, knurled handle.



CH06100

With lock.

CH06102

Without lock.



Titanium Delicate Touch® Micro Needle Holders

Straight. Round handle. Overall length 6in (15.2cm).



CH0393A

With lock.

CH0393

Without lock.

titanium



Titanium Delicate Touch® JACOBSON Micro Needle Holders

Round handle.



CH03184

Curved. With lock. Overall length 8in (20.3cm).

titanium



Delicate Touch® JACOBSON Micro Needle Holders

Round handle. Spiral cut. Mirror and Satin Finish. Overall length 8-1/4in (21cm).




◆ **CH06120**

Straight. With lock.

CH06122

Straight. Without lock.



CH06124

Curved. With lock.

CH06126

Curved. Without lock.



Delicate Touch® JACOBSON Micro Needle Holders

Flat handle. Mirror and Satin Finish.



CH06108

Straight. With lock. Overall length 8-1/2in (21.6cm).

CH06110

Straight. Without lock. Overall length 8-1/2in (21.6cm).

◆ **CH06158**

Straight. With lock. Overall length 9-1/2in (24.1cm).

CH06159

Straight. Without lock. Overall length 9-1/2in (24.1cm).



CH06112

Curved. With lock. Overall length 8-1/2in (21.6cm).

CH06114

Curved. Without lock. Overall length 8-1/2in (21.6cm).

CH06160

Curved. With lock. Overall length 9-1/2in (24.1cm).

CH06161

Curved. Without lock. Overall length 9-1/2in (24.1cm).



Delicate Touch® JACOBSON Micro Needle Holders

Short, round handle. Mirror and Satin Finish.
Overall length 9in (22.9cm).




CH06128

Straight. With lock.

CH06130

Straight. Without lock.



CH06132

Curved. With lock.

CH06134

Curved. Without lock.



Delicate Touch® CASTROVIEJO Needle Holders

Dust jaw surface for use with 6-0 to 8-0 suture. With lock.

CH8185

Straight.
Overall length
10in (25.4cm).



CH8187

Straight.
Overall length
12in (30.5cm).

CH8186

Curved.
Overall length
10in (25.4cm).

CH8188

Curved.
Overall length
12in (30.5cm).



Delicate Touch® JACOBSON Needle Holders

Dust jaw surface for use with 7-0 suture or smaller.
Round handle.

CH8210

With lock. Straight.
Overall length
10in (25.4cm).

CH8212

Without lock. Straight.
Overall length
10in (25.4cm).



CH8218

With lock. Straight.
Overall length
12in (30.5cm).

CH8220

Without lock. Straight.
Overall length
12in (30.5cm).

CH8214

With lock. Curved.
Overall length
10in (25.4cm).

CH8216

Without lock. Curved.
Overall length
10in (25.4cm).



CH8225

With lock. Curved.
Overall length
12in (30.5cm).

CH8230

Without lock. Curved.
Overall length
12in (30.5cm).



MILLS Microvascular Needle Holder

"Captain Hook."
Round BARRAQUER
handles. Delicate,
smooth,
curved jaws.



CH8752

Overall length
8in (20.3cm).



Vital® MILLS Vascular Needle Holders

Round BARRAQUER
handles. Spring handle.
Flat lock. Smooth
Tungsten Carbide
inserts. Overall length
8in (20.3cm).



CH8756

Straight.



CH8766

Curved.



Needle Holders

Porte-aiguilles
Nadelhalter
Porta-ago
Sostenes de agujas
Porta-agulhas
針ホルダー

Titanium Delicate Touch® JACOBSON Needle Holders

Finely tapered jaws. Round handle.

CH03170-S

With lock.
Straight.
Overall length
7-1/4in (18.4cm).



CH03172-S

Without lock.
Straight.
Overall length
7-1/4in (18.4cm)

CH03174-S

With lock.
Curved.
Overall length
7-1/4in (18.4cm).

CH03180-S

With lock.
Straight.
Overall length
8-1/4in (21cm).

CH03184-S

With lock.
Curved.
Overall length
8-1/4in (21cm).

CH03182-S

Without lock.
Straight.
Overall length
8-1/4in (21cm).

titanium



Titanium Delicate Touch® JACOBSON Micro Needle Holders

Round handle.

CH03170

Straight.
With lock.
Overall length
7-1/4in (18.4cm).



CH03174

Curved.
With lock.
Overall length
7-1/4in (18.4cm).

CH03172

Straight.
Without lock.
Overall length
7-1/4in (18.4cm).

CH03176

Curved.
Without lock.
Overall length
7-1/4in (18.4cm).

◆ CH03180

Straight.
With lock.
Overall length
8in (20.3cm).

CH03184

Curved.
With lock.
Overall length
8in (20.3cm).

CH03182

Straight.
Without lock.
Overall length
8in (20.3cm).



Titanium Delicate Touch® JACOBSON Micro Needle Holders

Flat handle.



CH03224

Straight.
With lock.
Overall length
7-1/4in (18.4cm).



CH03228

Curved.
With lock.
Overall length
7-1/4in (18.4cm).

CH03226

Straight.
Without lock.
Overall length
7-1/4in (18.4cm).

CH03230

Curved.
Without lock.
Overall length
7-1/4in (18.4cm).

CH03234

Straight.
With lock.
Overall length
8-1/2in (21.6cm).

CH03236

Curved.
Without lock.
Overall length
8-1/2in (21.6cm).

titanium



CH03224-S

With lock.
Straight.
Overall length
7-1/4in (18.4cm).

CH03226-S

Without lock.
Straight.
Overall length
7-1/4in (18.4cm).

CH03234-S

With lock.
Straight.
Overall length
8-1/2in (21.6cm).

CH03236-S

Without lock.
Straight.
Overall length
8-1/2in (21.6cm).

titanium



Titanium Delicate Touch® CASTROVIEJO Needle Holders

Dust jaw surface for 8-0 or smaller suture. With lock.



CH8180

Straight.
Overall length
9-1/2in (24.1cm).



CH8182

Curved.
Overall length
10in (25.4cm)

CH8181

Straight.
Overall length
10in (25.4cm)

CH8184

Curved.
Overall length
12in 30.5cm)

CH8183

Straight.
Overall length
12in (30.5cm)

titanium



Titanium Delicate Touch® JACOBSON Micro Needle Holders

Finely tapered jaws. Round handle.



CH0320012S

Straight jaws. With lock.
Overall length
12in (30.5cm).



CH0320412S

Curved jaws. With lock.
Overall length
12in (30.5cm).

CH03200-S

Straight jaws. With lock.
Overall length
10in (25.4cm).

CH03204-S

Curved jaws. With lock.
Overall length
10in (25.4cm).

CH03202-S

Straight jaws. Without lock.
Overall length
10in (25.4cm).

CH0320612S

Curved jaws. Without lock.
Overall length
12in (30.5cm).

titanium

CH0320212S

Straight jaws. Without lock.
Overall length
12in (30.5cm).




FINOCHIETTO Needle Holder

Straight jaws.



CH2460

Overall length
8in (20.3cm).



SAROT Needle Holders

Indented shanks. Straight jaws.




CH2415

Overall length
7-1/8in (18.1cm).



CH2410

Overall length
10-1/4in (26cm).



MASSON Needle Holder



CH2440

Overall length
10-1/2in (26.7cm).



Needle Holders

Porte-aiguilles
Nadelhalter
Porta-ago
Sostenes de agujas
Porta-agulhas
針ホルダー

Delicate Touch® Micro Ring-Handled Needle Holders



CH0620

Straight jaws.
Overall length
5in (12.7cm).

CH0624

Straight jaws.
Overall length
5-3/4in (14.6cm).

CH0628

Straight jaws.
Overall length
6-1/4in (15.9cm).

CH0632

Straight jaws.
Overall length
7-1/4in (18.4cm).

CH0636

Straight jaws.
Overall length
8-1/2in (21.6cm).

CH0640

Straight jaws.
Overall length
10in (25cm).

CH0634

Curved jaws.
Overall length
7-1/4" (18.4cm)

CH0638

Curved jaws.
Overall length
8-3/4in (22.2cm)

CH0642

Curved jaws.
Overall length
10in (25cm).



Titanium Delicate Touch® RICHARDS Micro Needle Holders

Extra-fine dust jaws
for use with 5-0 suture
and smaller. Straight.

CH03450

Overall length
6in (15.2cm).

titanium



Titanium Delicate Touch® RYDER Needle Holders

Dust jaw surface for use with
5-0 suture and smaller.
Straight.

CH03430

Overall length
6in (15.2cm).




CH03432

Overall length
6-3/4in (17.1cm).

CH03434

Overall length
7-1/2in (19cm).

titanium



Delicate Touch® Micro Needle Holders

RYDER-style.

CH0629

Overall length
6-1/4in (15.9cm).



CH0633

Overall length
7-1/4in (18.4cm).



Titanium Delicate Touch® MAYO-HEGAR Needle Holders

Dust jaw surface for use with
3-0 suture and smaller.
Straight.

CH03402

Overall length
6-1/4in (15.9cm).

CH03404

Overall length
7-1/8in (18.1cm).

CH03406

Overall length
8in (20.3cm).

CH03408

Overall length
8-5/8in (21.9cm).

titanium



Titanium Delicate Touch® Microvascular Needle Holders

Dust jaw surface for use with
5-0 suture and smaller.
Straight.



CH03422

Overall length
7-1/8in (18.1cm).

CH03424

Overall length
8in (20.3cm).

CH03426

Overall length
8-5/8in (21.9cm).

titanium



Titanium Delicate Touch® DeBAKEY Needle Holders

Dust jaw surface for use with
5-0 suture and smaller.
Straight.

CH03412

Overall length
7-1/8in (18.1cm).

CH03414

Overall length
8in (20.3cm).

CH03416

Overall length
8-5/8in (21.9cm).



CH03418

Overall length
10in (25.4cm).

titanium



Titanium Delicate Touch® Cardiovascular Needle Holders

Dust jaw surface for use with
5-0 suture and smaller.
RYDER-type jaws.
Straight.

CH03436

Overall length
7-1/8in (18.1cm).

CH03438

Overall length
8in (20.3cm).

CH03440

Overall length
8-5/8in (21.9cm).

CH03442

Overall length
10in (25.4cm).

titanium



Delicate Touch® JACOBSON Ring-Handled Micro Needle Holders

With lock.
Mirror and Satin Finish.



CH06148

Overall length
7-1/4in (18.4cm).

CH06150

Overall length
8in (20.3cm).

CH06152

Overall length
9in (22.9cm).



Forceps

Forceps
Zangen
Pinze
Fórcex
Fórcex
鉗子

JACOBSON Dressing Forceps

Thin handle, smooth jaws.



CH8515

Overall length
4in (10.2cm).



JACOBSON Dressing Forceps

Delicate. Broad handle.
Smooth jaws.



CH8520

Overall length
4-3/4in (12.1cm).



JACOBSON Dressing Forceps

Angled tips with tying platform.
Smooth jaws.



CH8525

Overall length
5-1/4in (13.3cm).



JACOBSON Dressing Forceps

Thin handle with stabilizing pin.
Smooth jaws.



CH8530

Overall length
6in (15.2cm).



JACOBSON Dressing Forceps

Thin, bayonet-shaped handles.
Smooth jaws.



CH8535

Overall length
6-1/2in (16.5cm).



JACOBSON Dressing Forceps

Thin handle with stabilizing pin.
Smooth jaws.



CH8540

Overall length
9in (22.9cm).



POTTS-SMITH Dressing Forceps

Fine, serrated tips.

CH5172

 Overall length
 7in (17.8cm).

CH5171

 Overall length
 8-1/4in (21cm).

CH5170

 Overall length
 9-1/2in (24.1cm).

CH5173

 Overall length
 10-1/4in (26cm).

CH5174

 Overall length
 12in (30.5cm).

CH5174-014

 Overall length
 14in (35.6cm).

POTTS-SMITH Tissue Forceps

Fine, serrated tips. 1 x 2 teeth.

CH5182

 Overall length
 7in (17.8cm).

CH5181

 Overall length
 8-1/4in (21cm).

CH5180

 Overall length
 9-1/2in (24.1cm).


CH5183

 Overall length
 10-1/4in (26cm).

CH5184

 Overall length
 12in (30.5cm).

CH5184-014

 Overall length
 14in (35.6cm).

Vital® POTTS-SMITH Forceps

 Delicate. Tungsten Carbide inserts.
 3,600 teeth per square inch.

CH5176

 Overall length
 7in (17.8cm).

CH5177

 Overall length
 8in (20.3cm).

CH5175

 Overall length
 9in (22.9cm).

CH5178

 Overall length
 9-1/2in (24.1cm).

CH5179

 Overall length
 10-1/4in (26cm).

KANTROWITZ Dressing Forceps

Curved, serrated tips. Angled jaws.


CH5195

 Overall length
 9-3/4in (24.7cm).


Forceps

Forceps
Zangen
Pinze
Fórcps
Fórcps
鉗子

GERALD Forceps

Straight. Flat handle.

CH0476

Serrated tips, 1mm wide.
Overall length
7in (17.8cm)

CH0477

Cross-serrated tips, 1mm wide.
Overall length
7-3/4in (19.7cm).



CH0479

Cross-serrated tips, 1mm wide.
Overall length
10in (25.4cm).

CH0479-12

Cross-serrated tips, 1.5mm wide.
Overall length
12in (30.5cm).



Titanium Delicate Touch® GERALD Forceps.

Straight. Flat handle.
1mm platform ring tips.



CH03750

Overall length
7in (17.8cm).

CH03752

Overall length
8in (20.3cm).

titanium



Titanium Delicate Touch® GERALD Forceps

Straight. Flat Handle.
Platform tips.



CH03760

Overall length
7in (17.8cm).

CH03762

Overall length
8in (20.3cm).

titanium



Delicate Touch® Micro Forceps

Lightweight. Straight.
1mm platform ring tips.
Round handle.



CH04282

Overall length
8-1/4in (21cm).



GERALD-DeBAKEY Forceps

Straight. Flat handle. Tips 1mm wide.




CH04242

Overall length
8in (20cm).



Delicate Touch® GERALD-DeBAKEY Forceps

Straight. Flat handle.
1mm wide DeBAKEY tips.



CH03740

Overall length
7in (17.8cm).

CH03742

Overall length
8in (20.3cm).



NELSON Lung Tissue Forceps

CH1500

6 x 7 teeth.
Overall length
9in (22.9cm).



RESANO-Type Forceps



CH5903

Overall length
7in (17.8cm).

CH5907

Overall length
9in (22.9cm).

CH5905

Overall length
10-1/4in (26.7cm).



DeBAKEY Dissecting Forceps

Atraumatic. Tips 1mm wide.
Round handle.

CH13178

Straight.
Overall length
7in (18cm).



V. Mueller Instrument Warranty

We warranty every medical device and accessory bearing the V. Mueller brand name to be free of functional defects in workmanship and materials when used normally for its intended use. Any V. Mueller product proving to be so defective will be replaced or repaired at no charge. Repairs or modifications performed by unauthorized personnel may void all product warranties and could affect performance and efficacy of the device.

Forceps

Forceps
Zangen
Pinze
Fórcps
Fórcps
鉗子

DeBAKEY-DIETRICH Vascular Tissue Forceps

Extra-delicate tips, 1mm wide.



◆ **CH5883**
Overall length
6in (15.2cm).

◆ **CH5885**
Overall length
7-3/4in (19.7cm).

CH5887
Overall length
9-1/2in (24.1cm).



DeBAKEY Vascular Tissue Forceps

Delicate jaws. Tips 1.5mm wide.

◆ **CH5893**
Overall length
6in (15.2cm).

◆ **CH5895**
Overall length
7-3/4in (19.7cm).


CH5897
Overall length
9-1/2in (24.1cm).



◆ **CH5899**
Overall length
12in (30.5cm).



DeBAKEY Vascular Tissue Forceps - Standard



◆ **CH5900**
Overall length
6-1/4in (15.9cm).


◆ **CH5902**
Overall length
7-3/4in (19.7cm).

◆ **CH5904**
Overall length
9-1/2in (24.1cm).



DeBAKEY Vascular Tissue Forceps

Heavy jaws taper from 2.5mm wide.



CH5900-1
Overall length
6-1/4in (15.9cm).

◆ **CH5902-1**
Overall length
7-3/4in (19.7cm).

CH5904-1
Overall length
9-1/2in (24.1cm).




Titanium DeBAKEY Tissue Forceps

CH5900-T
Tips 2mm wide.
Overall length
6in (15.2cm).

CH5902-T
Tips 2mm wide.
Overall length
7-3/4in (19.7cm).

CH5904-T
Tips 2mm wide.
Overall length
9-1/2in (24.1cm).

CH5906-T
Tips 2mm wide.
Overall length
12in (30.5cm).



CH5883-T
Tips 1mm wide.
Overall length
6in (15.2cm).

CH5885-T
Tips 1mm wide.
Overall length
7-3/4in (19.7cm).

CH5887-T
Tips 1mm wide.
Overall length
9-1/2in (24.1cm).

CH5912-T
Tips 3mm wide.
Overall length
6in (15.2cm).

CH5914-T
Tips 3mm wide.
Overall length
7-3/4in (19.7cm).


CH5916-T
Tips 3mm wide.
Overall length
9-1/2in (24.1cm).

CH5918-T
Tips 3mm wide.
Overall length
12in (30.5cm).

titanium



DeBAKEY Tissue Forceps



◆ CH5906

Jaws 2mm wide.
Overall length
12in (30.5cm).

CH5908

Jaws 2mm wide.
Overall length
14in (35.6cm).

CH5910


Jaws 2mm wide.
Overall length
16in (40.6cm).

CH5912

Heavy.
Jaws 2mm wide.
Overall length
18in (45.7cm).

CH5908-1

Jaws 3mm wide.
Overall length
14in (35.6cm).



CH5910-1

Jaws 3mm wide.
Overall length
16in (40.6cm).

CH5912-1

Heavy.
Jaws 3mm wide.
Overall length
18in (45.7cm).



DeBAKEY Tissue Forceps

Straight. Fluted handle.
Overall length
7-3/4in (19.7cm).



CH0444

1mm tips.

CH0446

2mm tips.



Titanium DeBAKEY Micro Forceps

Straight.
Round handle.
Tips 1mm wide.

CH03730

Overall length
7-1/8in (18.1cm).



CH03730-8

Overall length
8-1/4in (21cm).

titanium



DeBAKEY Vascular Tissue Forceps

Angled.
Jaws taper
from 2mm.

CH5920

Overall length
6in (15.2cm).

CH5922

Overall length
7-3/4in (19.7cm).

CH5924

Overall length
9-1/2in (24.1cm).


CH5926

Overall length
10in (25.4cm).



DeBAKEY Micro Tissue Forceps

Lightweight.
Round handle.
Tips 1mm wide.



CH04230

Overall length
7-1/8in (18.1cm).

CH04232

Overall length
8-1/4in (21cm).



Vital® DeBAKEY Tissue Forceps

Tips 2mm wide.



CH5200

Overall length
12in (30.5cm).



Forceps

Forceps
Zangen
Pinze
Fórcex
Fórcex
鉗子

Vital® DeBAKEY Needle Pulling Tissue Forceps

Flat handle. DeBAKEY jaws with Tungsten Carbide inserts.



CH5934

Tips 1.5mm wide.
Overall length
6in (15.2cm).

CH5930

Tips 2mm wide.
Overall length
6in (15.2cm).

CH5936

Tips 1.5mm wide.
Overall length
9-1/2in (24.1cm).

CH5935

Tips 1.5mm wide.
Overall length
7-3/4in (19.7cm).

CH5932

Tips 2mm wide.
Overall length
9-1/2in (24.1cm).

CH5931

Tips 2mm wide.
Overall length
7-3/4in (19.7cm).



Vital® DeBAKEY Micro Needle Pulling Tissue Forceps

Lightweight.
Round handle.
1.5mm wide
carbon inlay tips.



CH04266

Overall length
7-1/4in (18.4cm).



COOLEY Vascular Forceps

Delicate jaws.




CH6950-009

Overall length
5-1/2in (14cm).



COOLEY Vascular Tissue Forceps

COOLEY jaws.



CH5950

Jaws 1.5mm wide.
Overall length
6in (15.3cm).

CH5952

Jaws 2mm wide.
Overall length
7-3/4in (19.7cm).

CH5954

Jaws 2.3mm wide.
Overall length
9-1/2in (24.1cm).



COOLEY Vascular Tissue Forceps

1mm wide.




CH5940

Overall length
8in (20.3cm).



Cardiovascular Tissue Forceps

V. Mueller delicate jaws.


CH6000

 Overall length
 6-1/8in (15.6cm).

CH6001

 Overall length
 7-3/4in (19.8cm).

CH6002

 Overall length
 9-1/2in (24.1cm).

SELMAN Tissue Forceps

Fine, serrated, round tips.


CH6020

 Overall length
 6in (15.2cm).

CH6021

 Overall length
 8in (20.3cm).

CH6022

 Overall length
 9-3/4in (24.8cm).

WHEELER Forceps

CH6054


 Overall length
 6in (15.2cm).

CH6055

 Overall length
 7-1/8in (18.1cm).

CH6056

 Overall length
 8-3/4in (22.2cm).

ADAIR Tissue Forceps

CH1565

 Overall length
 9in (22.9cm).

ALLIS Tissue Forceps

 ◆ **CH1560**

 5 x 6 teeth.
 Overall length
 10in (25.4cm).

CH1560-MM

Mueller-Mirror® finish.

 ◆ **CH1562**

 5 x 6 teeth.
 Overall length
 12in (30.5cm).

CH1564

 5 x 6 teeth.
 Overall length
 14in (35.6cm).

CH1566

 5 x 6 teeth.
 Overall length
 16in (40.6cm).


CH1568

 5 x 6 teeth.
 Double-action
 box locks.
 Overall length
 18in (45.7cm).

BABCOCK Tissue Holding Forceps
CH1571

 Jaws 9.5mm wide.
 Overall length
 9-1/2in (24.1cm).

CH1570

 Jaws 16mm wide.
 Overall length
 9-1/2in (24.1cm).


Forceps

Forceps
Zangen
Pinze
Fórcps
Fórcps
鉗子

MAYO-RUSSIAN Tissue Forceps

Rugged, non-perforating teeth surround the perforated jaws. Deeply serrated handles.



CH1520

Overall length
14in (35.6cm).



Microsurgical Forceps

Counterbalanced. Round handle. Overall length 5-7/8in (15cm).



CH13159

Straight. Tip diameter 0.5mm.



Delicate Touch® Micro Forceps

Round handle. Straight. 1mm ring tips.

CH04216

Overall length
6in (15.2cm).



CH04220

Overall length
7-1/4in (18.4cm).

CH04224

Overall length
8-1/4in (21cm).

CH04226

Overall length
9in (22.9cm).


CH04228

Overall length
10in (25.4cm).



Titanium Delicate Touch® Micro Forceps

Straight. Round handle.



CH03158

Platform tips. Overall length
7-1/4in (18.4cm).



CH03160

1mm ring tips. Overall length
7-1/4in (18.4cm).

CH03168

1mm ring tips. Overall length
8in (20.3cm).

titanium



Delicate Touch® DENNIS Micro Forceps

Round handle. 1mm ring tips. Curved.



CH04221

Overall length
7-1/4in (18.5cm).



CH04225

Overall length
8-1/4in (21cm).



Delicate Touch® DENNIS Micro Forceps

Round handle. 2mm ring tips. Straight.



CH04320

Overall length
7in (17.8cm).



Titanium Delicate Touch® DENNIS Micro Forceps

Round handle.
1mm ring tips. Straight.

◆ CH03720

Overall length
7in (17.8cm).

CH03720-8

Overall length
8-1/4in (21cm).

CH03720-9

Overall length
9in (22.9cm).



CH03720-10

Overall length
10in (25.4cm).

titanium



Delicate Touch® Micro Forceps

Lightweight. 1 x 2 teeth.
Round handle.
Platform tips, 0.8mm wide.
Mirror and Satin Finish.

CH04235

Overall length
7-1/8in (18.1cm).



Delicate Touch® Micro Forceps

Lightweight. Straight. Platform tips.
Round handle. Mirror and Satin Finish.



CH04260

Overall length
7-1/8in (18.1cm).

CH04262

Overall length
8-1/4in (21cm).



MILLS Cardiovascular Tissue Forceps

Straight, 0.5mm ring tips.



◆ CH8770

Overall length
7in (17.8cm).



MILLS Cardiovascular Tissue Forceps

Round BARRAQUER
handles. Tear drop
ring tips, 0.5mm.



CH8770-001

Overall length
7in (17.8cm).



MILLS Thoracic Tissue Forceps

Round BARRAQUER
handles. Tear drop
ring tips, 1.5mm.



CH8772

Overall length
7in (17.8cm).



Forceps

Forceps
Zangen
Pinze
Fórcps
Fórcps
鉗子

Delicate Touch® Micro Forceps

Round handle.



CH0495

Straight.
Platform tips, 0.9mm wide.
Overall length
7-1/8in (18.1cm).

CH0497

Lightweight. Straight.
Overall length
7-1/8in (18.1cm).

CH0498

Straight.
Platform tips, 0.9mm wide.
Overall length
8-1/4in (21cm).

CH04204

Straight.
Platform tips, 0.9mm wide.
Overall length
10in (25.4cm).



CH0496

Curved.
Platform tips, 0.8mm wide.
Overall length
7-1/8in (18.1cm).

CH0497A

Lightweight. Curved.
Overall length
7-1/8in (18.1cm).

CH0499

Curved.
Platform tips, 0.8mm wide.
Overall length
8-1/4in (21cm).

CH04206

Curved.
Platform tips, 0.8mm wide.
Overall length
10in (25.4cm).



Delicate Touch® GERALD Forceps

Straight. Flat handle.



CH04246

1mm tips.
Overall length
7in (17.8cm).

CH04346

1mm platform tips.
Overall length
7in (17.8cm).

CH04244

1mm ring tips.
Overall length
7in (17.8cm).

CH04248

1mm tips.
Overall length
8in (20.3cm).

CH04249

1mm tips.
Overall length
9in (22.9cm).



Delicate Touch® Micro Forceps

Lightweight. Straight. Round handle. Satin and Mirror finish.



CH04250

1mm ring tips.
Overall length
7-1/8in (18.1cm).

CH04252

1mm ring tips.
Overall length
8-1/4in (21cm).

CH04254

2.5mm ring tips.
Overall length
7-1/8in (18.1cm).

CH04256

2.5mm ring tips.
Overall length
8-1/4in (21cm).



Delicate Touch® Micro Forceps

Lightweight.

Straight.

1mm platform ring tips.

Round handle.



CH04280

Overall length
7-1/8in (18.1cm).

CH04282

Overall length
8-1/4in (21cm).



LERICHE Hemostatic Forceps

Serrated jaws.
Overall length
6in (15.2cm).



CH7000

Delicate.
Straight.



CH7002

Straight.
1 x 2 teeth.




JACOBSON Mosquito Forceps

Very delicate.
Overall length
5in (12.7cm).



◆ **CH8608**

Straight jaws.



◆ **CH8610**

Curved jaws.




JACOBSON Hemostatic Forceps

Very delicate. Serrated jaws.
Overall length 7-1/2in (19.1cm).



◆ **CH8620**

Slightly curved. Overall length
7-1/2in (19.0cm).



CH8622

Fully curved. Overall length
7-1/4in (18.4cm).



CRILE-Type Clamps

Delicate. Serrated jaws,
length 1-1/2in (3.8cm).
Overall length
7-1/2in (19.1cm).




CH1705

Straight.



CH1707

Curved.




Forceps/Clamps

Forceps/Clamps
Zangen/Klemmen
Pinze/Morse
Fórceps/Pinças
Fórceps/Pinças
鉗子/クランプ


VANDERBILT Deep Vessel Forceps

Slender, curved,
serrated jaws.




◆ **CH1687**

Jaw length
1-3/4in (4.4cm).
Overall length
8-1/2in (21.6cm).



CRAFOORD-COLLIER Clamps

Curved. Delicate.



CH1690

Overall length
7-1/4in (18.4cm).

CH1694

Overall length
9-1/2in (24.1cm).




JULIAN Thoracic Artery Forceps

Slightly curved,
serrated jaws,
length 24mm.



CH1720

Overall length
9-1/4in (23.5cm).




RUMEL Thoracic Artery Forceps




CH1715

Slightly curved jaws.
Overall length
9-1/2in (24.1cm).




CH1716

Curved jaws.
Overall length
9-1/4in (23.5cm).



CH1717

Fully curved jaws.
Overall length
9-1/4in (23.5cm).




CH1718

Angled jaws.
Overall length
9-1/4in (23.5cm).



WOODWARD Thoracic Artery Forceps

Slightly curved jaws, length 32mm.




CH1719

Overall length
8-3/4in (22.2cm).



CRAFOORD Artery Forceps

Curved jaws.



CH1628

Overall length
9-1/2in (24.1cm).




SAROT Intrathoracic Artery Forceps

Serrated.




CH1680

Straight jaws,
length 3in (7.6cm).
Overall length
9-1/2in (24.1cm).



CH1682

Curved jaws,
length 3in (7.6cm).
Overall length
9-1/2in (24.1cm).



◆ **CH1685**

Curved jaws,
length 2-1/4in (5.7cm).
Overall length
9-3/4in (24.8cm).



Forceps

Forceps
Zangen
Pinze
Fórcps
Fórcps
鉗子

BRIDGE Deep Surgery Forceps

Serrated. Overall length 11in (27.9cm).




CH1724

Straight jaws,
length 39mm.




CH1724-001

Curved jaws,
length 36.5mm.



◆ CH1724-003

Delicate.
Curved jaws,
length 32mm.




V. Mueller Instrument Warranty

We warranty every medical device and accessory bearing the V. Mueller brand name to be free of functional defects in workmanship and materials when used normally for its intended use. Any V. Mueller product proving to be so defective will be replaced or repaired at no charge. Repairs or modifications performed by unauthorized personnel may void all product warranties and could affect performance and efficacy of the device.


BRIDGE Deep Surgery Forceps

Serrated. Curved. Double-action box locks.




CH1728

Overall length
18in (45.7cm).



HARRINGTON Forceps

Longitudinal serrations and
cross-serrated tip.



CH1740

Slightly curved jaws,
length 2-3/4in (7cm).
Overall length
12-1/2in (31.8cm).


CH1760

Right angled jaws,
length 2-3/8in (6cm).
Overall length
12in (30.5cm).



HARRINGTON-MIXTER Forceps

Right-angled jaws with longitudinal serrations and cross-serrated tip.




CH1770

Overall length
12-1/2in (31.8cm).




LAWRENCE Forceps

Curved.



CH1735

Overall length
11in (27.8cm).



KANTROWITZ Thoracic Clamps

Long handles.
Slender and delicate.



CH1721-001

Overall length
6in (15.2cm).

◆ **CH1721**
Overall length
7-1/2in (19.1cm).

CH1722
Overall length
9-1/2in (24.1cm).

CH1723
Overall length
10-3/4in (27.3cm).



KANTROWITZ Dressing Forceps

Curved, serrated tips. Angled jaws.

CH5197

Overall length
12in (30.5cm).



CH5199

Overall length
14in (35.6cm).




Forceps

Forceps
Zangen
Pinze
Fórcexs
Fórcexs
鉗子


MIXTER Thoracic Forceps (GEMINI-Type)

Delicate, right-angled, serrated jaws.



◆ CH1730-001

Overall length
7in (17.8cm).




CH1730-002

Overall length
8in (20.3cm).

◆ CH1730-003

Overall length
9in (22.9cm).



◆ CH1730-004

Overall length
11in (27.9cm).

CH1730-006

Overall length
14in (35.6cm).

CH1730-008

Overall length
16in (41cm).



MIXTER Forceps

Right-angled jaws with longitudinal serrations and cross-serrated tip.

◆ CH1725



Overall length
8-1/2in (21.6cm).

CH1726

Overall length
10-1/2in (26.7cm).


◆ CH1726-001

Overall length
11-1/2in (29.2cm).



UNIVERSITY OF MICHIGAN MIXTER Forceps

Right-angled, serrated jaws




CH1772

Overall length
9-1/4in (23.5cm).



FINOCHIETTO Thoracic Forceps

Longitudinally serrated right-angled jaws.
Perforated tip.




CH1780

Overall length
9-3/4in (24.8cm).



DeBAKEY-COOLEY Cardiovascular Forceps

Longitudinal serrations with cross-serrated tips.
Slightly curved jaws.




CH1775

Overall length
10in (25.4cm).



PENNINGTON Tissue Holding Forceps

Triangular serrated jaws.




CH1590

Overall length
8in (20.3cm).



LOVELACE Lung Forceps

Angled shanks.
Jaws 1in (25mm) wide.



CH1610

Overall length
7-3/8in (18.7cm).




Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ


SAROT Bronchus Clamps

Jaws 1-1/2in (38mm) wide.
Overall length 9in (22.9cm).




◆ **CH1650**
Curved right.

CH1660
Curved left.



LEES Bronchus Clamp

Jaw length 1-1/2in (4cm).




CH1675

Overall length
9-1/4in (23.5cm).




LEES Right Angle Clamp

Jaw length 1-1/2in (4cm).



CH1670

Overall length
8in (20.3cm).



Customer support

V. Mueller products are sold only by authorized V. Mueller representatives. All instruments and equipment of our manufacture bear the V. Mueller trademark.

For more information, please contact your sales representative, send an email to **GMB-VM-Europe@carefusion.com** or visit **carefusion.com/vmueller**.

Due to changing surgical trends, some products described in this catalog may be discontinued or not available in your country.

Your V. Mueller sales representative may be able to suggest a suitable substitute if one is available.


DeBAKEY-SEMB Ligature Carrier

Fully curved.


CH1809Overall length
10-1/4in (25.6cm).**SEMB Dissecting Forceps/Ligature Carriers**

Perforated tips for use as ligature carrier.


Overall length 9-1/2in (24.1cm).

**CH1800**

Slightly curved jaws.

**CH1805**

Fully curved jaw.

**FAVALORO-SEMB Ligature Carriers****CH1806**Half-curved jaws.
Overall length
11-5/16in (28.7cm).**CH1808**Fully curved jaws.
Overall length
11-1/8in (28.2cm).

Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ



Approximator Clamps

Without frame. Double.

Artery	Vein	Size	Color
CH12000	CH12002	0.4mm-1mm	Black
CH12001	CH12003	0.4mm-1mm	Black
CH12008	CH12010	0.6mm-1.5mm	Black
CH12009	CH12011	0.6mm-1.5mm	Black
CH12016	CH12018	1mm-2.25mm	Black
CH12017	CH12019	1mm-2.25mm	Black



Approximator Clamps

With frame. Double.

Artery	Vein	Size	Color
CH12004	CH12006	0.4mm-1mm	Black
CH12005		0.4mm-1mm	Black
CH12012	CH12014	0.6mm-1.5mm	Black
CH12013	CH12015	0.6mm-1.5mm	Black
CH12020	CH12022	1mm-2.25mm	Black
CH12021	CH12023	1mm-2.25mm	Black



Single Clamps

Without frame. Packaged in pairs.

Artery	Vein	Size	Color
CH11000	CH11002	0.4mm-1mm	Black
CH11001	CH11003	0.4mm-1mm	Black
CH11004	CH11006	0.6mm-1.5mm	Black
CH11005	CH11007	0.6mm-1.5mm	Black
CH11008	CH11010	1mm-2.25mm	Black
CH11009	CH11011	1mm-2.25mm	Black



Double Clamps

With frame. Stationary.

Artery	Vein	Size
CH11012		0.4mm-1mm
CH11016		0.6mm-1.5mm
CH11020	CH11022	1mm-2.25mm

General Purpose Clamps


Packaged in pairs.

CH13029Curved.
Vessel size 0.3mm–1mm.
Black.**CH13030**Curved.
Vessel size 0.6mm–1.4mm.**CH13155**Angled.
Vessel size 1mm–2.25mm.**Hand-Applied Clamps**Flat jaws. Single. Without frame.
Packaged in pairs.**CH13036**

Vessel size 2mm–5mm.

CH13037Vessel size 2mm–5mm.
Black.**CH13063**

Vessel size 1.5mm–3.5mm.

CH13064Vessel size 1.5mm–3.5mm.
Black.**Hand-Applied Approximator Clamps**Flat jaws. Double.
Without frame.**CH13034**Vessel size
2mm–5mm.**CH13035**Vessel size
2mm–5mm.
Black.**CH13061**Vessel size
1.5mm–3.5mm.**CH13062**Vessel size
1.5mm–3.5mm.
Black.**Clamp Applying Forceps**

Without lock. Angled.

CH13144Overall length
5-7/8in (15cm).**Customer support**

V. Mueller products are sold only by authorized V. Mueller representatives. All instruments and equipment of our manufacture bear the V. Mueller trade-mark.

For more information, please contact your sales representative, send an email to GMB-VM-Europe@carefusion.com or visit carefusion.com/vmueller.

Due to changing surgical trends, some products described in this catalog may be discontinued or not available in your country.

Your V. Mueller sales representative may be able to suggest a suitable substitute if one is available.

Bulldog Clamps

Clamps bulldog
Bulldogklemmen
Serracavo
Pinzas de mordaza
Pinça Bulldog
ブルドッグクランプ

WELDON Miniature Bulldog Clamps

Finely serrated jaws,
length 1/2in (1.3cm).
Overall length
1-3/8in (3.5cm).



◆ **CH6077-001**

Straight jaws.

◆ **CH6077-002**

Curved jaws.



Micro Cross-Action Bulldog Clamps

Finely serrated jaws.



◆ **CH6079**

Angled jaws,
length 1/2in (1.3cm).
Overall length
2-1/8in (5.4cm).

◆ **CH6078**

Angled jaw,
length 3/8in (1cm).
Overall length
2in (5.2cm).

CH6078-001

Straight jaws,
length 3/8in (1cm).
Overall length
2in (5.2cm).

◆ **CH6079-001**

Straight jaws,
length 1/2in (1.3cm).
Overall length
2-1/4in (5.7cm).



DIETHRICH Angled Bulldog Clamp

Serrated.



CH6076

Jaw 20mm.
Overall length
1-7/8in (4.8cm).



Bulldog Clamps

CH6143

Angled jaws,
length 3/8in (1cm).
Overall length
1-1/4in (3.3cm).

CH6144


Angled jaws,
length 1/2in (1.3cm).
Overall length
1-1/2in (3.8cm).

CH6145

Angled jaws,
length 3/4in (1.8cm).
Overall length
1-1/2in (3.8cm).

CH6148

Straight jaws,
length 3/4in (1.8cm).
Overall length
1-1/2in (3.8cm).



COOLEY Bulldog Clamps

Adjustable screw to regulate jaw tension.
Straight, finely serrated jaws.

CH6030

Jaw length
3/4in (1.9cm).
Overall length
2in (5.1cm).

CH6031

Jaw length
1-1/8in (2.9cm).
Overall length
2-5/8in (6.7cm).

CH6032

Jaw length
1-7/8in (4.8cm).
Overall length
3-3/8in (8.6cm).


CH6033

Jaw length
2-3/8in (6cm).
Overall length
5in (12.7cm).



COOLEY Bulldog Clamps

Adjustable screw to regulate jaw tension.
Curved, finely serrated jaws.



◆ CH6035

Jaw length 3/4in (1.9cm).
Overall length
1-3/4in (4.4cm).

CH6036

Jaw length 1in (2.5cm).
Overall length
2-1/2in (6.4cm).

CH6037

Jaw length 1-1/2in (3.8cm).
Overall length
3-1/8in (7.9cm).

CH6038

Jaw length 2-3/8in (6cm).
Overall length
4-1/2in (11.7cm).



COOLEY Cross-Action Bulldog Clamps


Curved jaw.

CH6040

Jaw length 5/8in (1.6cm).
Overall length
2-7/8in (7.3cm).

CH6041

Jaw length 1in (2.5cm).
Overall length
3-1/4in (8.3cm).



CH6042

Jaw length 1-1/2in (3.8cm).
Overall length
3-7/8in (9.8cm).

CH6043

Jaw length 2-1/4in (5.7cm).
Overall length
4-1/2in (11.4cm).



COOLEY Cross-Action Bulldog Clamps

Straight jaw.

CH6045

Jaw length
7/8in (2.2cm).
Overall length
3in (7.6cm).



CH6047

Jaw length
1-7/8in (4.8cm).
Overall length
4-1/4in (10.8cm).

CH6048

Jaw length
2-5/8in (6.7cm).
Overall length
5in (12.7cm).

CH6046

Jaw length
1-1/8in (2.9cm).
Overall length
3-3/8in (8.6cm).



CH6050

Straight jaw,
length 3cm.
Overall length
2-1/2in (6.5cm).

CH6060

Straight jaw,
length 4.8cm.
Overall length
3-1/8in (8cm).

CH6052

Curved jaw,
length 3cm.
Overall length
2-1/2in (6.5cm).

CH6062

Curved jaw,
length 4.8cm.
Overall length
3-1/8in (8cm).



V. Mueller® Bulldog Clamps

V. Mueller delicate jaws.

Bulldog Clamps

Clamps bulldog
Bulldogklemmen
Serracavo
Pinzas de mordaza
Pinça Bulldog
ブルドッグクランプ

GREGORY Bulldog Clamps

1 x 2 teeth.

CH6089

Straight jaw,
length 1-3/16in (3cm).
Overall length
3-1/2in (8.9cm).

CH6091

Straight jaw,
length 1-3/4in (4.4cm).
Overall length
4-1/4in (10.8cm).

CH6094

Curved jaw,
length 1-15/16in (4.9cm).
Overall length
3-1/2in (8.9cm).


CH6090

Curved jaw,
length 1/2in (1.25cm).
Overall length
4-1/4in (10.8cm).



GREGORY Soft Bulldog Clamps

Overall length 4-3/4in (12cm).



CH6093-001


Curved right.

CH6093-002

Curved left.



GREGORY DeBAKEY Bulldog Clamp




CH6092

Angled 45°.
Jaw length 7cm.
Overall length
4-3/4in (12.1cm).



Cross-Action Bulldog Clamps

V. Mueller delicate jaws, straight.



CH6100

Length 1.7cm.
Overall length
3in (7.6cm).

CH6101

Length 3cm.
Overall length
3-3/8in (8.5cm).


CH6102

Length 4.8cm.
Overall length
4-1/8in (10.5cm).



Cross-Action Bulldog Clamps

V. Mueller delicate jaws. Curved.



CH6110

Length 2.3cm.
Overall length
3in (7.5cm).

CH6111

Length 3.2cm.
Overall length
3-3/8in (8.5cm).

CH6113

Length 7cm.
Overall length
4-3/4in (12cm).



Customer support

V. Mueller products are sold only by authorized V. Mueller representatives. All instruments and equipment of our manufacture bear the V. Mueller trademark.


For more information, please contact your sales representative, send an email to GMB-VM-Europe@carefusion.com or visit carefusion.com/vmueller.

Due to changing surgical trends, some products described in this catalog may be discontinued or not available in your country.

Your V. Mueller sales representative may be able to suggest a suitable substitute if one is available.

JACOBSON Bulldog Clamps

Adjustable screw to control tension.
Finely serrated jaws, length 7/8in (2.1cm).
Overall length 1-3/4in (4.5cm).



◆ CH8505

Straight.

◆ CH8510

Curved.



GLOVER Bulldog Clamps

Straight. DeBAKEY jaws.
Adjusting screw to regulate tension.



CH6065

Jaw length 1.8cm.
Overall length
2in (5.1cm).

CH6066

Jaw length 2.7cm.
Overall length
2-1/2in (6.4cm).

CH6067

Jaw length 4cm.
Overall length
3-1/2in (8.9cm).

CH6068

Jaw length 6cm.
Overall length
4-3/4in (12.1cm).




Bulldog Clamps

Clamps bulldog
Bulldogklemmen
Serracavo
Pinzas de mordaza
Pinça Bulldog
ブルドッグクランプ

GLOVER Bulldog Clamps

Curved. DeBAKEY jaws.
Adjusting screw to regulate tension.



CH6070

Jaw length 1.8cm.
Overall length
2in (5.1cm).

CH6071

Jaw length 2.7cm.
Overall length
2-1/2in (6.4cm).

CH6072

Jaw length 4cm.
Overall 3-1/2in (8.9cm).


CH6073

Jaw length 6cm.
Overall length
4-3/4in (12.1cm).



GLOVER Bulldog Clamp

COOLEY jaws. Angled 90°.



CH6075

Jaw length 3.1cm.
Overall length
1-5/8in (4.1cm).



DeBAKEY Cross-Action Bulldog Clamps

Straight jaw.



CH6081

Jaw length 2cm.
Overall length
3in (7.6cm).

CH6082

Jaw length 3cm.
Overall length
3-3/8in (8.6cm).


CH6083

Jaw length 4.5cm.
Overall length
4-1/8in (10.5cm).



DeBAKEY Cross-Action Bulldog Clamps

Curved jaw.



CH6085

Jaw length 2cm.
Overall length
3in (7.6cm).

CH6086

Jaw length 3cm.
Overall length
3-3/8in (8.7cm).

CH6087

Jaw length 4.5cm.
Overall length
4-1/8in (10.5cm).

CH6088


Jaw length 6.5cm.
Overall length
4-3/4in (12cm).



Clamps bulldog/Forceps
Bulldogklemmen/Zangen
Serracavo/Pinze
Pinzas de mordaza/Fórceps
Pinça Bulldog/Fórceps
ブルドッグクランプ/鉗子

Titanium Bulldog Clamps

Lightweight, corrosion-resistant.



Strong Close

Straight	Curved	Jaw Length
CH6116	CH6118	3cm
CH6122	CH6124	3.5cm
CH6128	CH6130	4cm
	CH6140	5cm

Slight Close

Straight	Curved	Angled	Jaw Length
CH6115	CH6117	CH6120	3cm
CH6121	CH6123	CH6126	3.5cm
	CH6129	CH6132	4cm

titanium



Bulldog Clamp Applying Forceps

CH5160

Overall length
9-1/4in (23.5cm).



DeBAKEY Ring-Handled Bulldog Clamps

Overall length
5-1/8in (13cm).




CH6095

Straight.



CH6096

Curved jaws,
length 4.8cm.




Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

DeBAKEY Ring-Handled Bulldog Clamp

Overall length
5in (12.7cm).




DeBAKEY Ring-Handled Bulldog Clamp

Overall length
3-3/4in (10cm).






DeBAKEY Clamp

Curved shanks.






DeBAKEY Miniature Multi-Purpose Clamps

Jaw length 1-1/8in (2.9cm). Overall length 6-1/2in (16.5cm).





SELMAN Peripheral Blood Vessel Forceps


Straight COOLEY jaws.

**CH6150**Jaw length 1.3cm.
Overall length
5in (12.7cm).**CH6152**Jaw length 2.5cm.
Overall length
5-1/2in (14cm).**SELMAN Peripheral Blood Vessel Forceps**


Angled COOLEY jaws.

**CH6151**Jaw length 1.3cm.
Overall length
5in (12.7cm).**CH6153**Jaw length 2.1cm.
Overall length
5-1/2in (14cm).**SELMAN Peripheral Blood Vessel Forceps**

Fully curved COOLEY jaws.

**CH6170**Jaws
15mm x 8.5mm deep.
Overall length
5in (12.7cm).**COOLEY Vascular Clamp**

Extra small. Curved shanks.


Angled jaws,
length 1-1/2in (3.8cm).◆ **CH6950-003**Overall length
4-7/8in (12.4cm).

Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

COOLEY Vascular Clamp

Extra small. Straight shanks.
Fully curved jaws,
length 1in (2.5cm).




CH6950-004

Overall length
4-1/2in (11.4cm).



COOLEY Vascular Clamp

Extra small. Straight shanks.
Partial occlusion jaws,
length 1-1/8in (2.9cm).




◆ **CH6950-005**

Overall length
4-5/8in (11.7cm).



COOLEY Vascular Clamp

Extra small. Straight shanks.
Straight jaws,
length 1-1/2in (3.8cm).




CH6950-001

Overall length
5in (12.7cm).



COOLEY Vascular Clamp

Extra small. Curved shanks.
Angled jaws,
length 1-3/4in (4.4cm).




◆ **CH6950-006**

Overall length
5in (12.7cm).



COOLEY Vascular Clamp

Extra small.
Angled shanks.
Straight jaws,
length 1-1/2in (3.8cm).

**CH6950-002**

Overall length
5in (12.7cm).

**COOLEY Clamps**

Jaw length 3cm.

CH6516

Straight.
Overall length
5-1/4in (13.3cm).

CH6517


Angled jaws.
Overall length
5-1/4in (13.3cm).

CH6518


Straight jaws.
Angled shanks.
Overall length
5-1/4in (13.3cm).

CH6519

Curved jaws.
Overall length
5-1/2in (14cm).

**COOLEY-SATINSKY Clamp**

Jaw length 30mm,
depth 5mm.

**CH6521**

Overall length
5-1/2in (14cm).

**COOLEY Patent Ductus Clamps**

Jaw length
1-1/4in (3.2cm).
Overall length
6-1/4in (15.9cm).

**CH6180**

Straight shanks.

CH6181

Angled shanks.




Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

COOLEY Anastomosis Clamps

Double-angled jaws are calibrated in 5mm increments.
Angled shanks. Overall length 6-1/2in (16.5cm).



CH6567

Jaw length
3/4in (1.9cm),
depth 7.3mm.


CH6568

Jaw length
3/4in (1.9cm),
depth 6mm.



COOLEY Coarctation Clamp

Angled shanks.
Jaw length 1-1/2in (3.8cm).




CH6183

Overall length
6-3/4in (17.1cm).




COOLEY Coarctation Clamps

Straight shanks. Jaw length 1-1/2in (3.8cm).
Overall length 6-3/4in (17.1cm).




CH6182

Straight jaw.



CH6185

Curved jaw.



Customer support

V. Mueller products are sold only by authorized V. Mueller representatives. All instruments and equipment of our manufacture bear the V. Mueller trademark.


For more information, please contact your sales representative, send an email to GMB-VM-Europe@carefusion.com or visit carefusion.com/vmueller.

Due to changing surgical trends, some products described in this catalog may be discontinued or not available in your country.


Your V. Mueller sales representative may be able to suggest a suitable substitute if one is available.

COOLEY Tangential Occlusion Clamps


Double-angled jaws are calibrated in 5mm increments.
Overall length 7-3/4in (19.7cm).

**CH6570**

Jaw length 3/4in (1.9cm),
depth 6mm.

**CH6572**


Jaw length 1in (2.5cm),
depth 7mm.

**CH6574**

Jaw length 1-3/8in (3.5cm),
depth 10mm.

**SATINSKY Clamp**

COOLEY jaws.
Jaw length 1in (2.5cm),
depth 5mm.

**CH6520**


Overall length
5-1/2in (14cm).

**KITZMILLER-Type Clamps**

COOLEY jaws.
Overall length
5-1/4in (13.4cm).

CH6510

Angled left.

**CH6512**

Angled right.



Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

Customer support

V. Mueller products are sold only by authorized V. Mueller representatives. All instruments and equipment of our manufacture bear the V. Mueller trademark.

For more information, please contact your sales representative, send an email to **GMB-VM-Europe@carefusion.com** or visit **carefusion.com/vmueller**.

Due to changing surgical trends, some products described in this catalog may be discontinued or not available in your country.

Your V. Mueller sales representative may be able to suggest a suitable substitute if one is available.

MULLER Clamps

V. Mueller delicate jaws.
Jaw length 2.1cm.
Overall length
7-1/4in (18.4cm).



CH6500

Straight shanks.


CH6502

Angled shanks.



CALMAN Closure Clamp

Angled handles.
Long jaw, 1mm thick.




CH7902

Overall length
4-3/8in (11.1cm).



THOMPSON Closure Clamp

Handles slightly curved.
Long jaw, 1mm thick.




CH7910

Overall length
6-1/8in (15.5cm).



KARTCHNER Internal Clamp


Long, slightly curved handles.
Long jaw, length 1-3/8in (3.5cm).

**CH7912**

Overall length
6-3/8in (16.2cm).

**COOLEY Peripheral Vascular Clamp**


Light. Straight. Jaw length 1-3/4in (4.4cm).

**CH6190**

Overall length
6-3/4in (17.1cm).

**COOLEY Peripheral Vascular Clamps**

Light. Jaws angled 55°.

**CH6198**

Jaw length 1-3/8in (3.5cm).
Overall length
6-1/2in (16.5cm).

CH6199

Jaw length 1-3/4in (4.4cm).
Overall length
6-3/4in (17.1cm).

CH6201


Jaw length 1-3/4in (4.4cm).
Overall length
7in (17.8cm).

CH6200

Jaw length 2-3/4in (7cm).
Overall length
7in (17.8cm).

**COOLEY Patent Ductus Clamps**

Straight jaws,
length 1-1/8in (3cm).
Overall length
8-1/4in (21cm).

**CH6320**

Straight shanks.

CH6325

Angled shanks.




Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

LELAND JONES Peripheral Vascular Clamps

Jaw length 3-1/2in (8.9cm). Overall length 8-1/4in (21cm).



CH7190

Straight.

CH7192

Jaws angled 15°.

CH7194


Jaws angled 30°.



DeBAKEY Patent Ductus and Peripheral Vascular Clamps

Straight jaws, length 1-1/8in (3cm).

Overall length 8in (20.3cm).



CH7100

Straight shanks.

CH7102


Angled shanks.



DeBAKEY Curved Peripheral Vascular Clamp


A small, lightweight clamp.

Curved jaws,
length 1-7/8in (4.8cm).




CH7155

Overall length
6in (15.2cm).




DeBAKEY Peripheral Vascular Clamp

Lightweight shanks. Slight curve in shank with full angle
in jaw distal to box lock. Jaw length 2-3/8in (6cm).




◆ **CH7140**

Overall length
7in (17.8cm).



DeBAKEY-BAINBRIGE Vascular Clamps**CH7116**


Curved. Jaw length 1-1/2in (4cm).
Overall length 6-7/8in (17.4cm).

**CH7117**


Curved. Jaw length 2-3/8in (6cm).
Overall length 7-3/4in (19.7cm).

**DeBAKEY S-Shaped Peripheral Vascular Clamp**

Curved jaws and light, springy handles conform to an S-shape.
Jaw length 3-3/8in (8.5cm).

**CH7170**


Overall length
8in (20.3cm).

**GREGORY Profunda Clamps**◆ **CH7221**

Medium.
Jaw length
2-1/2in (6.2cm).
Overall length
7in (18cm).

CH7222


Large.
Jaw length
2-7/8in (7.3cm).
Overall length
7-3/4in (19.7cm).

◆ **CH7220**

Small.
Jaw length
1-3/4in (4.4cm).
Overall length
5-7/8in (13.7cm).

COOLEY-BAUMGARTEN Aortic Clamp

Slightly curved jaws and shanks.
Jaw length 2-1/2in (6.4cm).

**CH6215**

Overall length
7-1/2in (19.1cm).




Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

BAILEY Aortic Clamps

Angled shafts.



CH7284


Extra small jaws, length 1-1/8in (29mm), depth 8mm. Overall length 8in (20.3cm).

CH7286

Small jaws length 1-3/4in (44mm), depth 11mm. Overall length 8in (20.3cm).

CH7288


Small jaws length 2-1/8in (54mm), depth 11mm. Overall length 9-3/4in (24.7cm).



HOWARD Aortic Aneurysm Clamps

Jaw length 10.5cm

Overall length 11-3/4in (29.8cm).




CH7206

Angled.
DeBAKEY 2 x 3 jaw.

CH7202

Reverse curve shanks.
DeBAKEY 2 x 3 jaw.




CH7204

Angled.
DeBAKEY 1 x 2 jaw.



ZANGER Abdominal Aorta Clamps


Angled.

CH7214Overall length
9-1/2in (24.1cm).**CH7216**Overall length
12in (30.5cm).


Instrument sizes, particularly overall lengths, are specified in this catalog as accurately as possible. Slight variations may occur due to metal crafting and finishing operations. Refer to the introductory pages for clarification on measurements. We reserve the right to modify products and/or services depicted herein without customer notice.

ZANGER Abdominal Aorta Clamps

Acutely-curved shanks.

CH7210Overall length
9-1/2in (24.1cm).**CH7212**Overall length
12-1/4in (31cm).**DeBAKEY Aortic Clamp**

Curved jaws and light, springy shanks conform to S-shape.


**CH7175**Jaw length 3in (7.6cm).
Overall length
7-1/4in (18.4cm).

Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ


DeBAKEY Aorta Clamps

Slightly curved jaws.



CH7245

Jaw length
3-3/4in (9.5cm).
Overall length
10-1/2in (26.5cm).




◆ **CH7247**

Jaw length 4in
(10cm).
Overall length
12in (30.5cm).



MORRIS Aorta Clamps

Jaw length 2-1/2in (6.4cm) from angle.



CH7510

Overall length
8-1/8in (20.6cm).





CH7512

Overall length
9-1/2in (24.1cm).



FITZGERALD Aorta Clamps

Slightly curved COOLEY jaws.

CH6260Jaw length
3-1/4in (8.5cm).
Overall length
10in (25.5cm).**CH6262**Jaw length
3-3/4in (9.5cm).
Overall length
10-3/4in (26.5cm).**CH6270**Jaw length
4-1/2in (11.2cm).
Overall length
12-1/4in (31cm).**CH6272**Jaw length
5-1/4in (13.2cm).
Overall length
13in (33cm).**WYLIE Abdominal Aorta Clamp****CH7456**Overall length
8-1/2in (21.6cm).**DIETHRICH Aorta Clamp**1 x 2 vascular teeth.
Jaw length 1-3/4in (45mm),
depth 5/8in (15mm).**CH7280**Overall length
8-3/8in (21.3cm).**LEMOLE-STRONG Aorta Clamps**Jaw width 1-3/4in (4.4cm),
depth 7/8in (2.2cm).**CH7490**Straight.
Overall length
7-3/4in (19.6cm).**CH7492**Angled.
Overall length
8in (20.3cm).

Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

COOLEY Aorta Clamp

Curved shanks.




CH6442

Jaw length 2-7/8in (7.3cm), depth 3/4in (1.9cm). Overall length 10in (25.4cm).

COOLEY Aortic Aneurysm Clamp

Slightly curved jaws.




CH6280

Jaw length 3in (7.6cm). Overall length 10-3/4in (27.3cm).

COOLEY Aorta Clamps

Curved shanks.



CH6444

Jaw length 1-7/8in (4.7cm), depth 1/2in (1.3cm). Overall length 8in (20.3cm).


CH6446

Jaw length 2-1/8in (5.4cm), depth 1/2in (1.3cm). Overall length 10-3/4in (27.3cm).




KAY Clamps


DeBAKEY jaws. Overall length 8in (20.3cm).

**CH7350**

Jaw length 1-3/4in (4.5cm), depth 1-1/4in (3.2cm).

**CH7351**

Jaw length 1-3/4in (4.5cm), depth 3/4in (1.9cm).




Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ


PARSONNET Aorta Clamps

The jaws, available in four sizes. Overall length 10-1/4in (26cm).




CH7355-001

Jaw length at base
1-1/8in (2.8cm).
Overall length
9-3/8in (23.8cm).




CH7355-002

Jaw length at base
1-1/4in (3.2cm).
Overall length
9-3/8in (23.8cm).



CH7355-003

Jaw length at base
1-1/2in (3.8cm).
Overall length
9-1/2in (24.1cm).




CH7355-004

Jaw length at base
1-3/4in (4.4cm).
Overall length
9-3/4in (24.7cm).




DeBAKEY "Sidewinder" Aorta Clamps

Double curved jaws.




CH7252

Jaw length 1-3/4in (4.5cm).
Overall length
9in (22.9cm).




CH7251

Jaw length 2-3/8in (6cm).
Overall length
10-1/2in (26.7cm).






DERRA Miniature Anastomosis Clamps

Spring Handle.

CH6814Small. Jaw 1/2in (1.5cm).
Overall length
6-7/8in (17.4cm).**CH6816**Medium. Jaw 3/4in (2cm).
Overall length
7-1/8in (18cm).**CH6818**Large. Jaw 1in (2.5cm).
Overall length
6-3/4in (17.2cm).**COOLEY-DERRA Clamps**

DERRA-type jaws in 5mm increments. Overall length 6-1/2in (16.5cm).


**CH6233**Jaw length 3/4in (1.9cm),
depth 3/8in (1cm).**CH6234**Jaw length 7/8in (2.2cm),
depth 1/2in (1.3cm).**CH6235**Jaw length 1-1/4in (3.2cm),
depth 5/8in (1.6cm).

Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ


DeBAKEY SATINSKY Clamps

Overall length 10in (25.4cm).




CH7300

Small jaws,
length 2in (5.1cm),
depth 5/16in (0.9cm).



CH7302

Medium jaws, length 2-1/8in (5.4cm),
depth 3/4in (1.9cm).



CH7304

Large jaws, length 2-3/8in (6cm),
depth 3/4in (1.9cm).



DERRA Vestibular Clamp

Longitudinally serrated jaws.

CH5098

Jaws 1-3/4in x 3in
(44mm x 75mm) wide.
Overall length
9-1/2in (24.1cm).



Customer support

V. Mueller products are sold only by authorized V. Mueller representatives. All instruments and equipment of our manufacture bear the V. Mueller trademark.





For more information, please contact your sales representative, send an email to GMB-VM-Europe@carefusion.com or visit carefusion.com/vmueller.

Due to changing surgical trends, some products described in this catalog may be discontinued or not available in your country.

Your V. Mueller sales representative may be able to suggest a suitable substitute if one is available.

COOLEY Coarctation Clamps

Straight jaws.

CH6300Straight shanks.
Jaw length 2in (5.1cm).
Overall length
8-3/4in (22.2cm).**CH6302**Straight shanks.
Jaw length 2-3/4in (7cm).
Overall length
9-3/4in (24.8cm).◆ **CH6310**Angled shanks.
Jaw length 2in (5.1cm).
Overall length
8-3/4in (22.2cm).**CH6312**Angled shanks.
Jaw length 2-3/4in (7cm).
Overall length
9-3/4in (24.8cm).**CRAFOORD Coarctation Clamp**Longitudinally serrated
jaws, length 2-3/4in (7cm).**CH5120**Overall length
9-1/4in (23.5cm).**CRAFOORD Coarctation Forceps****CH6360**Slightly curved jaws,
length 2-3/4in (7cm).
Overall length
9-1/4in (23.5cm).**DeBAKEY Coarctation Vascular Clamps**Straight, non-crushing jaws, length 2-1/8in (5.4cm).
Overall length 9in (22.9cm).**CH7110**

Straight shank.

CH7112

Angled shanks.




Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

GLOVER Coarctation Clamps

CH7180

Extra-long. Straight.
Jaw length 2-3/4in (7cm).
Overall length 10-1/2in (26.7cm).



CH7182

Extra-long. Angled.
Jaw length 2-3/4in (7cm).
Overall length 10-1/2in (26.7cm).

CH7184

Straight. Jaw length 2-1/8in (5.4cm).
Overall length 9in (23cm).


CH7186

Angled. Jaw length 2-1/8in (5.4cm).
Overall length 8-7/8in (22.5cm).




CASTANEDA Multi-Purpose Clamps

Jaw length 1-3/8in (3.5cm). Overall length 5in (12.7cm).




CH6803

Straight.



CH6800

Jaws angled 30°.




◆ CH6801

Jaws angled 50°.




CH6802

Jaws angled 90°.



**DeBAKEY-CASTANEDA
Clamps**

Delicate, narrow. 1.4mm DeBAKEY jaws.
Overall length 4-3/4in (12.1cm).

**CH7813**


Jaws angled 45°.
Curved shaft.

CH7814


Jaws angled 60°.
Curved shaft.

CH7815


Jaws angled 90°.
Curved shaft.

**DeBAKEY Multi-Purpose Clamps**


Jaws angled 90°.

**CH7262**


Jaw length 2in (5.1cm).
Overall length
8in (20.3cm).

**CH7264**


Jaw length 2-1/2in (6.4cm).
Overall length
8-1/2in (21.6cm).

**DeBAKEY Multi-Purpose Clamps**

Curved shanks. Jaws angled 60°.

**CH7272**

Jaw length
2in (5.1cm).
Overall length
8-1/4in (21cm).

**CH7274**

Jaw length
2-1/2in (6.4cm).
Overall length
8-1/2in (21.5cm).

CH7276

Jaw length
4in (10.2cm).
Overall length
12in (30.5cm).



Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

DeBAKEY Multi-Purpose Clamp

Jaws angled 25°, length 1-1/2in (3.8cm).




CH7268

Overall length
9in (22.9cm).

DeBAKEY Multi-Purpose Aorta Clamp

Jaw angled 90°. Jaw length 4-1/4in (10.8cm).



CH7265

Overall length
9in (22.9cm).

DeBAKEY Multi-Purpose Clamp

Curved shanks. Jaws angled 60°, length 1-1/2in (3.8cm).




CH7270

Overall length
9-1/4in (23.5cm).

COOLEY Multi-Purpose Clamps

Jaws in 5mm increments.



CH6210

Jaw angled 30°,
length 1-1/2in (3.8cm).
Overall length
6-1/2in (16.5cm).

CH6211


Jaw angled 60°,
length 1-1/8in (2.9cm).
Overall length
6-1/8in (15.6cm).

CH6212


Jaw angled 90°,
length 1-1/2in (3.8cm).
Overall length
5-1/2in (14cm).

COOLEY Multi-Purpose Curved Clamps

Acutely curved jaws in 5mm increments.
Angled shanks. Overall length 8in (20.3cm).

**CH6620**


Jaw length 1-1/4in (3.2cm),
depth 5/8in (1.6cm).

**CH6622**


Jaw length 1-1/2in (3.8cm),
depth 3/4in (1.9cm).

**COOLEY Multi-Purpose Angled Clamps**


Jaws angled 60° half the distance from the distal tip.

**CH6660**


Jaw length 1-1/2in (3.8cm).
Overall length 9-1/2in (24.1cm).

**CH6661**

Jaw length 2-1/8in (5.4cm).
Overall length 9in (22.9cm).

**CH6662**

Jaw length 2-3/4in (7cm).
Overall length 9-1/2in (24.1cm).




Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ


CASTANEDA Anastomosis Clamps

Overall length 6in (15.2cm).




CH6807

Jaw length 5/8in (1.5cm), depth 5mm.



CH6808

Jaw length 3/4in (1.8cm), depth 6mm.




CH6809

Jaw length 7/8in (2.1cm), depth 7mm.



COOLEY Double-Angled Clamps

Jaws in 5mm increments, angled 90°.




CH6230

Jaw length 1 in (2.5cm), depth 6.8mm.
Overall length 6-1/2in (16.5cm).




COOLEY Double-Angled Clamps

Jaws in 5mm increments, angled 60°.




CH6225

Jaw length 1-1/2in (3.8cm), depth 7.6mm.
Overall length 7in (17.8cm).




CH6226

Jaw length 1-1/4in (3.2cm), depth 7.6mm.
Overall length 7-1/4in (18.4cm).





COOLEY Carotid, Subclavian, and Renal (CSR) Artery Clamps

Jaws angled 1/3 of the distance from the distal tip and in 5mm increments.

**HENLY Artery Clamp**

Double-angled jaws.

**DeBAKEY Artery Clamp**Ratchet handle. Curved with 1 x 9 ratchets.
Jaw length 2-1/2in (6.5cm).

Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

COOLEY Partial Occlusion Clamps

CH6650

Jaw length
1in (2.5cm),
depth 7mm.
Overall length
8in (20.3cm).

CH6652

Jaw length
2in (5.1cm),
depth 1.2cm.
Overall length
9in (22.9cm).

CH6654


Jaw length
2-7/8in (7.3cm),
depth 1.2cm.
Overall length
10-3/4in (27.3cm).

CH6651

Jaw length
1-3/8in (3.5cm),
depth 7mm.
Overall length
8-1/2in (21.6cm).


CH6653

Jaw length
2-3/8in (6cm),
depth 1.2cm.
Overall length
9-1/2in (24.1cm).



V. Mueller® Auricular Appendage Clamp

Jaw length 4.9cm, base 2.4cm, depth 1.3cm.




CH6422

Overall length
10in (25.4cm).



COOLEY Auricular Appendage Forceps

Jaws 1-5/8in (4.1cm) x 3/8in (1cm).




CH6648

Overall length
10in (25.4cm).



COOLEY Curved Clamp

Light. Fully curved jaws.



**CH6220**Jaw length 1-1/8in (2.9cm).
Overall length
6-1/4in (15.9cm).**Customer support**

V. Mueller products are sold only by authorized V. Mueller representatives. All instruments and equipment of our manufacture bear the V. Mueller trademark.

For more information, please contact your sales representative, send an email to **GMB-VM-Europe@carefusion.com** or visit **carefusion.com/vmueller**.

Due to changing surgical trends, some products described in this catalog may be discontinued or not available in your country.

Your V. Mueller sales representative may be able to suggest a suitable substitute if one is available.

DeBAKEY Tangential Occlusion Clamps◆ **CH7290**Jaw length 3.5cm,
depth 1.1cm.
Overall length
7-3/4in (19.7cm).**CH7294**Jaw length 5.5cm,
depth 1.5cm.
Overall length
10-1/2in (26.7cm).**CH7298**Jaw length 7.5cm,
depth 1.7cm.
Overall length
11in (28cm).**CH7292**Jaw length 4.8cm,
depth 1.6cm.
Overall length
10-1/4in (26cm).**CH7296**Jaw length 6cm,
depth 1.4cm.
Overall length
10-3/4in (27.3cm).

Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

COOLEY Iliac Clamps

Curved shanks, angled 45°.
Straight jaws, length 2-3/8in (6cm).



CH6640

Overall length
8-1/4in (21.0cm).

CH6642

Overall length
10in (25.4cm).




V. Mueller Instrument Warranty

We warranty every medical device and accessory bearing the V. Mueller brand name to be free of functional defects in workmanship and materials when used normally for its intended use. Any V. Mueller product proving to be so defective will be replaced or repaired at no charge. Repairs or modifications performed by unauthorized personnel may void all product warranties and could affect performance and efficacy of the device.


COOLEY Curved Clamps

Acutely curved jaws are in 5mm increments.



CH6236

Jaw length 7/8in (2.2cm),
depth 3/8in (1cm).
Overall length
6in (15.2cm).




CH6237

Jaw length 1-3/8in (3.5cm),
depth 1/2in (1.3cm).
Overall length
6-1/4in (15.9cm).






COOLEY Curved Cardiovascular Clamp

Curved jaws, length 2-1/8in (5.4cm).

**CH6355**Overall length
8-1/4in (21cm).**COOLEY Cardiovascular Clamp**


Curved jaws. Jaw length 2in (5.1cm).

**CH6646**Overall length
10-3/4in (27.3cm).**DeBAKEY Curved Clamp**Strongly curved jaws,
length 2-3/8in (6cm).**CH7230**Overall length
10-1/4in (26cm).**DeBAKEY Curved Clamp**Acutely curved jaws,
length 3-3/8in (8.6cm).**CH7240**Overall length
10-1/2in (26.7cm).

Clamps


Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

DeBAKEY-BAHNSON Curved Clamps



CH7231

Curved jaws,
length 2-1/2in (6.4cm).
Overall length
10-3/4in (27.3cm).




CH7232

Slightly curved jaws,
length 3in (7.6cm).
Overall length
11-1/4in (28.6cm).



DeBAKEY Hypogastric Clamp



CH7340

Overall length
8-1/2in (21.5cm).




WYLIE "J" Clamps

Acutely curved jaws,
length 1-1/4in (3.3cm), depth 1in (2.5cm).

CH7450

Overall length
8-1/2in (21.6cm).




CH7452


Overall length
9-5/8in (24.4cm).

CH7454

Overall length
9-1/8in (23.2cm).



WYLIE Hypogastric Clamp




CH7458

Overall length
9-1/2in (24.1cm).



DALE Clamps**CH7472**


Small. Jaw length 1/2in (1.5cm).
Overall length 6-7/8in (17.5cm).

**CH7474**

Medium.
Jaw length 3/4in (2cm).
Overall length 7in (17.8cm).

**DALE Improved Clamp**


Jaw length 1-1/2in (4cm).

**CH7470**

Overall length
7-1/2in (19cm).

**DEMOS Artery Clamp**


Jaw length 1-1/2in (4cm).

**CH7480**


Overall length
8-1/2in (21.5cm).

**CHERRY Clamp**

Handle angled 60°.
2 x 3 teeth.
Curved jaws,
length 5-3/4in
(14.6cm).

**CH7460**

Overall length
12-1/2in (31.7cm).




Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

HENDREN Cardiovascular Clamps

Very Delicate. COOLEY jaws in 5mm increments. Overall length 8-1/4in (21cm).



CH6540

Jaw base 6mm.

CH6542

Jaw base 12mm.

CH6544


Jaw base 18mm.

CH6546

Jaw base 24mm.


CH6548

Jaw base 30mm.



GLOVER Curved Forceps

Curved jaws, length 1-3/4in (4.5cm), depth 7mm.




CH7115

Overall length
8-1/4in (21cm).




GLOVER Right Angle Curved Clamp

Right-angled jaws, length 2in (5cm), depth 7mm.




CH7188

Overall length
7-3/4in (19.7cm).




GLOVER Spoon Anastomosis Clamp


Jaw length 2-7/16in (6.2cm).

**CH6340**Overall length
8-1/2in (21.6cm).**HENDREN Cardiovascular Clamps**

Delicate. COOLEY jaws in 5mm increments. Overall length 6-1/4in (15.9cm).

**CH6522**



Jaw base 12mm.

**CH6524**

Jaw base 18mm.

**GLOVER Curved Clamps**

Curved COOLEY jaws.

**CH6350**Jaw length 2in (5.1cm).
Overall length
8-1/2in (21.6cm).**CH6352**Jaw length 2-3/4in (7cm).
Overall length
9-3/4in (24.8cm).

Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

GLOVER Spoon-Shaped Forceps



CH7118

Curved jaw, length 1-5/8in (4cm), depth 7mm.
Overall length 8-1/2in (21.6cm).


CH7119

Curved jaw, length 1-7/8in (4.8cm), depth 8mm.
Overall length 9-1/2in (24.1cm).



GLOVER Auricular Appendage Clamp

Jaw length 1-3/4in (4.3cm), depth 7mm.




CH7310

Overall length 9in (22.9cm).




HARKEN Clamps

Curved jaw with cross-serrations.




CH6464

Jaw length 2-3/8in (6cm).
Overall length 9-3/8in (23.8cm).



CH6462

Jaw length 2-7/8in (7.4cm).
Overall length 9-1/2in (24.1cm).




CH6461

Jaw length 3-1/8in (8cm).
Overall length 9-3/4in (24.9cm).



HARKEN Clamps

Curved jaws with COOLEY teeth.
Overall length 10in (25.4cm).

**COOLEY Caval Occlusion Clamps**

Clamps encourage flow through indwelling, cannula, or catheter.
Jaw ring clamps around catheter. Angled shanks.

CH7970

Ring size 24Fr. Overall length
6-1/2in (16.5cm).

CH7972

Ring size 28Fr. Overall length
7-1/4in (18.4cm).

CH7974

Ring size 32Fr. Overall length
7-1/4in (18.4cm).

CH7974-001


Ring size 34Fr. Overall length
7-1/4in (18.4cm).

CH7975

Ring size 36Fr. Overall length
7-1/4in (18.4cm).

**COOLEY Bronchus Clamp**


Designed with COOLEY jaws.
Jaw length 2in (5.1cm).

**CH6400**


Overall length
8-1/4in (21cm).

**COOLEY Bronchus Clamp**

Angled shanks.
COOLEY jaws angled 90°,
length 3-3/4in (9.5cm).

**CH6405**

Overall length
9-1/2in (24.1cm).



Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

JAVID Carotid Artery Bypass Shunt Clamps

Designed to hold tubes in place during vascular procedures.



◆ CH7950

Angled jaws close to approximately 4mm diameter ring at distal tip.
Overall length
6-1/4in (15.9cm).



◆ CH7952

Angled jaws close to approximately 6mm diameter ring at distal tip.
Overall length
7in (17.8cm).



Fogarty Hydragrip® Clamps

For vascular use. Jaws accept 33mm inserts.

◆ CV5020

Straight. Overall length
6-1/2in (16.5cm).

◆ CV5025

Angled. Overall length
4-5/8in (11.8cm).



Fogarty® Spring Clip Applicator Forceps

◆ CV5110

For use with 6mm Fogarty® spring clips.
Overall length
7-1/4in (18.4cm).

CV5111

For use with 12mm Fogarty® spring clips.
Overall length
7in (17.8cm).



Fogarty Hydragrip® Clamps

Jaws accept 61mm inserts.

◆ CV5030

Straight. Overall length
11-1/2in (29.2cm).

◆ CV5040

Angled. Overall length
8-3/4in (22.2cm).



Fogarty Hydragrip® Clamp

Jaws accept 61mm inserts.

◆ **CV5045**Angled.
Overall length
8-1/2in (21.6cm).**Fogarty Hydragrip® Clamp**

Jaws accept 86mm inserts.

◆ **CV5050**Angled.
Overall length
9-1/4in (23.5cm).**Fogarty Hydragrip® Forceps**

Jaws accept 86mm inserts.

◆ **CV5055**Straight.
Overall length
12-1/4in (31.1cm).**COSGROVE Slim Fit™ Clamp**Has a more slender clamp head.
Inserts not included.◆ **CV1033**Fogarty-type jaw,
length 33mm.
Working length 31cm.◆ **CV1061**Fogarty-type jaw,
length 61mm.
Working length 33cm.

Clamps

Clamps
Klemmen
Morse
Pinzas
Pinças
クランプ

COSGROVE Quick Bend™ Flex Clamp

Has a shortened neck that bends closer to the jaw for greater flexibility. Inserts not included.

◆ **CV1133**

Fogarty-type jaws accept 33mm inserts.

◆ **CV1161**

Fogarty-type jaws accept 61mm inserts.

◆ **CV1186**

Fogarty-type jaws accept 86mm inserts.



DeBAKEY Vascular Tunneler

Tip diameter 16mm.
Mirror and Satin finish.

CH12900

Working length
16-1/2in (42cm).
Overall length
21-1/4in (54cm).



NOON (Modified) Vascular Access Tunneler

Includes three interchangeable olive tips: 6, 7 and 9mm diameters.
Mirror finish.


CH12920

Overall length
9-1/4in (23.5cm).



HALLMAN Tunneling Forceps

Double-action. Cross-serrated ring jaws.



CH1900

Overall length
14-5/8in (37.1cm).



Chest Tube Passers

◆ CH04185

Without ratchet.
Overall length
9-3/4in (24.8cm).



CH04186

With ratchet.
Overall length
9-3/4in (24.8cm).



CH04189

Without ratchet.
Overall length
14in (35.6cm).

V. Mueller® Blunt Hooks

Matte and Mirror finish.
Overall length 7in (17.8cm).



CH8672

2mm rounded tip.



CH8673

1.3mm ball tip.



Punches

Emporte-pièces
Stanzen
Punzoni
Sacabocados
Punções
穿孔

MERZ Aorta Punch

Straight. Punch size 5mm.



CH5625

Overall length
8-1/2in (21.6cm).



Customer support

V. Mueller products are sold only by authorized V. Mueller representatives. All instruments and equipment of our manufacture bear the V. Mueller trademark.


For more information, please contact your sales representative, send an email to **GMB-VM-Europe@carefusion.com** or visit **carefusion.com/vmueller**.

Due to changing surgical trends, some products described in this catalog may be discontinued or not available in your country.

Your V. Mueller sales representative may be able to suggest a suitable substitute if one is available.

NEW ORLEANS Endarterectomy Stripper Set

Five stripping loops, sizes 4mm, 5mm, 6mm, 7mm and 9mm on 12in (30.5cm) handle. Handle is threaded on proximal end. Additional 12in (30.5cm) female-threaded handle included to extend stripping loops.



CH8320



LEATHER (Modified MILLS) Retrograde Valvulotomes

Blade on inside edge of tip. Tip 2.9mm. Mirror finish.



CH8680

Overall length
9-1/2in (24.1cm).



CH8680-001


Overall length
13in (33cm).



Aiguille/Canules
Nadel/Kanülen
Ago/Cannule
Aguja/Cánulas
Agulha/Cânulas
針/カニューレ

LEATHER Antegrade Valvulotome

Cutting surface on leading edge. Mirror finish.

**CH8681**

Overall length
7in (17.8cm).

**ZANGER Cannula**

Malleable. Luer lock. Tip 4mm.

CH8084

Overall length
2-1/2in (6.3cm).

**Infusion Cannula**

Malleable. Tip 1.5mm.

**CH8085**

Overall length
2-1/8in (5.4cm).

**JACOBSON Vessel Knife**

Pointed, concave blade.

**CH8586**

Overall length
7-1/4in (18.4cm).



Aiguille/Canules
Nadel/Kanülen
Ago/Cannule
Aguja/Cánulas
Agulha/Cânulas
針/カニューレ

Heparin Needles

13GA. Luer lock.
Mirror finish.

CH8080

Overall length
2in (5.1cm).

CH8083

Overall length
3in (7.6cm).



CLAGETT "S" Cannulas

Double curved, Luer hub.
Mirror finish.

CH3390

13GA.
Overall length
5-1/2in (14cm).

CH3391

14GA.
Overall length
5-1/2in (14cm).


CH3392

14GA.
Overall length
3-1/8in (7.9cm).



PATTERSON Trocars

Mirror and Satin finish.




CH2806

Small.
Sheath 5.9mm O.D.
Overall length
7-1/2in (19cm).




CH2804

Medium.
Sheath 8.4mm O.D.
Overall length
7-1/2in (19cm).



CH2802

Medium-large.
Sheath 10mm O.D.
Overall length
8in (20.3cm).



CH2800

Large.
Sheath 11.4mm O.D.
Overall length
8in (20.3cm).



JACOBSON Endarterectomy Spatulas

Semi-sharp blade, 10mm x 1mm wide.
Overall length 7in (17.8cm).



CH8580

Straight.



CH8585

Angled 45°.



Vascular Suction Tubes

CH8040

Long suction basket tip.
Working length from
angle to distal tip
3-1/2in (8.9cm).
Overall length
11in (28.0cm).

CH8042

Small detachable tip.
Tip angled 135°.
Overall length
10in (25.4cm).



FELL Suction Tube

Tip 9mm O.D.

CH8050

Overall length
13in (33cm).



COOLEY Vascular Suction Tube

Round,
perforated end
7.5mm O.D.

CH8057

Working length
7in (17.8cm).
Overall length
12in (30.5cm).



COOLEY Graft Suction Tube

Open end 8mm O.D.,
small opening on side.

CH8055

Working length
8-1/2in (21.6cm).
Overall length
14in (35.6cm).



COOLEY Intracardiac Suction Tube

Round, perforated end
9mm O.D.

CH8059

Working length
7in (17.8cm).
Overall length
12in (30.5cm).



COOLEY Sump Suction Tubes

Tip diameter 7.5mm.
Overall length
13in (33cm).

CH8060

Angled 30°,
left ventricle.



Suction Tubes

Tubes d'aspiration
Absaugschläuche
Tubi di aspirazione
Tubos de succión
Tubos de sucção
吸引チューブ

Vascular Suction Tube

Working length from angle to distal tip
5in (12.7cm). Overall length
8-1/4in (21cm).

CH8044

Long, fluted tip.

CH8046

Tips tapers with 1.7mm hole.

CH8047

Tip tapers with 2mm hole.

CH8048

Tip tapers with 3mm hole.

CH8049

Tip tapers with 4mm hole.



Small Suction Tubes

No. 1, fluted, fixed tip.

CH3300

Working length
2-1/2in (6.4cm).
Overall length
7-1/2in (19.1cm).



Small Suction Tubes

No. 2, perforated, removable tip.

CH3302

Working length
2-1/2in (6.4cm).
Overall length
7-1/2in (19.1cm).



V. Mueller's VITAL® Advantage


Helping you choose the right instrument for the right procedure is a specialty of V. Mueller, the market leader in specialty surgical instruments.

Vital® Needle Holder User's Guide provides three easy ways to choose the most appropriate Needle Holder for your surgical needs.

1. **By suture size.** Identifying the appropriate needle holder to use with the suture and needle selected. This is easy with three separate charts – 4-0, 5-0-6-0 and 7-0 and finer.
2. **By alphabetical cross-reference** of instrument name to catalog number.
3. **By Needle Table** of sutures and suture needles offered by two leading manufacturers.

As a surgical specialist, we're sure you'll agree – starting with the right needle holder is one of your surest ways to reduce repairs and prolong the life of your Vital® instrument.

For 5-0 and 6-0 Suture Size




Cat. No.	Description	Length
10-20mm Needle Length – Needle Table B		
CH2411	SAROT	10-3/8" (26.3cm)
CH2416	SAROT	7-1/8" (18.1cm)
CH2417	ROCHESTER, delicate	7-1/8" (18.1cm)
CH2420	JULIAN, delicate, straight	8" (20.3cm)
CH2421	JULIAN, delicate, curved	8" (20.3cm)
CH2422	CRILE-WOOD	7-1/8" (18.1cm)
CH2423	CRILE-WOOD, curved	7-1/8" (18.1cm)
CH2424	CRILE-WOOD	8" (20.3cm)
CH2425	CRILE-WOOD	9" (22.9cm)
CH2426	CRILE-WOOD, curved	9-1/4" (23.5cm)
CH2426-001	CRILE-WOOD, left handed	9-1/4" (23.5cm)
CH2427	CRILE-WOOD	10-3/8" (26.3cm)
CH2428	CRILE-WOOD	11-3/4" (29.8cm)
CH2429	CRILE-WOOD	14in (35.6cm)

Needle Tables

Table A				Table B				Table C				
4-0 and Larger Suture				5-0 and 6-0 Suture				7-0 and Finer Suture				
.22-65mm Length				10-20mm Length				2-11mm Length				
Ethicon			Davis & Geck	Ethicon		Davis & Geck		Ethicon		Davis & Geck		
BB	G	SC	C	PRE	BB	MO	C	PE	BV	P	C	LE
BP	J	SH	CE	SBE	BV	MH	CE	PR	C	RB	CE	LO
C	KS	SHB	CS	SCE	C	P	CO	PRE	CC	SIF	CF	PR
CCS	LR	ST	CT	SLO	CC	PS	D	SBE	CIF	ST	CO	PRE
CP	LS	STP	CV	T	CCS	RB	DO	SCE	CS	TG	CV	SE
CPS	M	TE	D	TE	CPS	RD	DT	SE	CTC	VAS	CTE	T
CPX	MH	TS	DO	TH	CT	S	DTE	T	G	V	DT	TE
CT	MO	UCL	MT	TS	CTX	SH	LA	TE			DTO	TS
CTB	OS	UR	NE	TT	FN	SM	LE	Y			LA	V
CTX	P	V	PR		FS	TF	LO	YE			LB	YE
EN	PC	X			G	TG	MR					
FS	PS	XLH			M	V						
FSL	PSL											
FSLX	S											


Recommended VITAL® Needle Holders

For up to 4-0 Suture Size



Cat. No.	Description	Length
22-65mm Needle Length – Needle Table A		
CH2433	WANGENSTEEN	10-3/4" (27.3cm)
CH2435	MAYO-HEGAR	11-3/4" (29.8cm)
CH2437	MAYO-HEGAR	14in (35.6cm)
CH2438	MAYO-HEGAR	18in (45.7cm)
CH2442	MASSON	10-3/8" (26.3cm)
CH2461	FINOCHIETTO	8" (20.3cm)
CH2471	FINOCHIETTO	10-1/2" (26.7cm)

For 7-0 and Finer Suture Size



Cat. No.	Description	Length
2-11mm Needle Length – Needle Table C		
CH2406	COOLEY Microvascular	6" (15.2cm)
CH2407	COOLEY Microvascular	6-3/4" (17.1cm)
CH2408	COOLEY Microvascular	8" (20.3cm)
CH2419	CRILE-WOOD	6" (15.2cm)

code	page	class	code	page	class	code	page	class	code	page	class
CH03158	66	I	CH03432	56	I	CH04262	67	I	CH0628	56	I
CH03160	66	I	CH03434	56	I	CH04266	64	I	CH0629	56	I
CH03168	66	I	CH03436	57	I	CH04280	68	I	CH06310	50	I
CH03170	54	I	CH03438	57	I	CH04282	68	I	CH0632	56	I
CH03170-S	54	I	CH03440	57	I	CH04282	60	I	CH0633	56	I
CH03172	54	I	CH03442	57	I	CH04320	66	I	CH06330	50	I
CH03172-S	54	I	CH03450	56	I	CH04346	68	I	CH0634	56	I
CH03174	54	I	CH03720	67	I	CH0444	63	I	CH06340	50	I
CH03174-S	54	I	CH03720-10	67	I	CH0446	63	I	CH0636	56	I
CH03176	54	I	CH03720-8	67	I	CH0476	60	I	CH0638	56	I
CH03180	54	I	CH03720-9	67	I	CH0477	60	I	CH0640	56	I
CH03180-S	54	I	CH03730	63	I	CH0479	60	I	CH0642	56	I
CH03182	54	I	CH03730-8	63	I	CH0479-12	60	I	CH0650	50	I
CH03182-S	54	I	CH03740	61	I	CH0495	68	I	CH0652	50	I
CH03184	52	I	CH03742	61	I	CH0496	68	I	CH0654	50	I
CH03184-S	54	I	CH03750	60	I	CH0497	68	I	CH0656	50	I
CH03184-S	54	I	CH03752	60	I	CH0497A	68	I	CH0658	50	I
CH0320012S	55	I	CH03760	60	I	CH0498	68	I	CH0660	50	I
CH03200-S	55	I	CH03762	60	I	CH0499	68	I	CH0662	50	I
CH0320212S	55	I	CH0393	52	I	CH06100	52	I	CH0664	50	I
CH03202-S	55	I	CH0393A	52	I	CH06102	52	I	CH0670	51	I
CH0320412S	55	I	CH04185	125	I	CH06108	52	I	CH0672	51	I
CH03204-S	55	I	CH04186	125	I	CH06110	52	I	CH0674	51	I
CH0320612S	55	I	CH04189	125	I	CH06112	52	I	CH0676	51	I
CH03224	54	I	CH04204	68	I	CH06114	52	I	CH0682	51	I
CH03224-S	54	I	CH04206	68	I	CH06120	52	I	CH0684	51	I
CH03226	54	I	CH04216	66	I	CH06121	51	I	CH0686	51	I
CH03226-S	54	I	CH04220	66	I	CH06122	52	I	CH0688	51	I
CH03228	54	I	CH04221	66	I	CH06123	51	I	CH0690	51	I
CH03230	54	I	CH04224	66	I	CH06124	52	I	CH0692	51	I
CH03234	54	I	CH04225	66	I	CH06125	51	I	CH0694	51	I
CH03234-S	54	I	CH04226	66	I	CH06126	52	I	CH0696	51	I
CH03236	54	I	CH04228	66	I	CH06127	51	I	CH07274	33	I
CH03236-S	54	I	CH04230	63	I	CH06128	53	I	CH07276	33	I
CH03402	56	I	CH04232	63	I	CH06130	53	I	CH07420	32	I
CH03404	56	I	CH04235	67	I	CH06132	53	I	CH07422	32	I
CH03406	56	I	CH04242	61	I	CH06134	53	I	CH07424	32	I
CH03408	56	I	CH04244	68	I	CH06148	57	I	CH07426	32	I
CH03412	57	I	CH04246	68	I	CH06150	57	I	CH07428	32	I
CH03414	57	I	CH04248	68	I	CH06152	57	I	CH07450	32	I
CH03416	57	I	CH04249	68	I	CH06158	52	I	CH07452	32	I
CH03418	57	I	CH04250	68	I	CH06159	52	I	CH07454	32	I
CH03422	57	I	CH04252	68	I	CH06160	52	I	CH07456	32	I
CH03424	57	I	CH04254	68	I	CH06161	52	I	CH07458	32	I
CH03426	57	I	CH04256	68	I	CH0620	56	I	CH07460	33	I
CH03430	56	I	CH04260	67	I	CH0624	56	I	CH07460-1	33	I

code	page	class	code	page	class	code	page	class	code	page	class
CH07462	33	I	CH1130	30	Ila	CH1265	28	Ila	CH1610	75	I
CH07462-1	33	I	CH1131	30	Ila	CH1275-1	27	Ila	CH1628	71	Ila
CH07464	33	I	CH1132	30	Ila	CH1280-1	27	Ila	CH165	15	I
CH07464-1	33	I	CH1133	30	Ila	CH1285-1	27	Ila	CH1650	76	Ila
CH07466	33	I	CH1140-1	25	Ila	CH1287	25	Ila	CH1660	76	Ila
CH07466-1	33	I	CH120	15	I	CH1288	25	Ila	CH1670	76	Ila
CH07468	33	I	CH12000	78	Ila	CH1289	25	Ila	CH1675	76	Ila
CH07468-1	33	I	CH12001	78	Ila	CH1290	27	Ila	CH1680	71	I
CH07470	33	I	CH12002	78	Ila	CH12900	125	I	CH1682	71	I
CH07472	33	I	CH12003	78	Ila	CH1290-1	27	Ila	CH1685	71	I
CH0758	38	I	CH12004	78	Ila	CH12920	125	I	CH1687	70	I
CH0764	38	I	CH12005	78	Ila	CH130	15	I	CH1690	70	Ila
CH0812	41	I	CH12006	78	Ila	CH13029	79	Ila	CH1694	70	Ila
CH0814	41	I	CH12008	78	Ila	CH13030	79	Ila	CH1705	69	Ila
CH1000	26	I	CH12009	78	Ila	CH13034	79	Ila	CH1707	69	Ila
CH1005	26	I	CH12010	78	Ila	CH13035	79	Ila	CH1715	70	I
CH1008	26	I	CH12011	78	Ila	CH13036	79	Ila	CH1716	70	I
CH1009	26	I	CH12012	78	Ila	CH13037	79	Ila	CH1717	70	I
CH1010	26	I	CH12013	78	Ila	CH13061	79	Ila	CH1718	70	I
CH1030	30	Ila	CH12014	78	Ila	CH13062	79	Ila	CH1719	71	I
CH1030-001	30	Ila	CH12015	78	Ila	CH13063	79	Ila	CH1720	70	I
CH1030-002	30	Ila	CH12016	78	Ila	CH13064	79	Ila	CH1721	73	Ila
CH1030-003	30	Ila	CH12017	78	Ila	CH13144	79	I	CH1721-001	73	Ila
CH1030-004	30	Ila	CH12018	78	Ila	CH13155	79	Ila	CH1722	73	Ila
CH1030-005	30	Ila	CH12019	78	Ila	CH13159	66	I	CH1723	73	Ila
CH1030-006	30	Ila	CH12020	78	Ila	CH13178	61	I	CH1724	72	I
CH1030-007	30	Ila	CH12021	78	Ila	CH1340	29	Ila	CH1724-001	72	I
CH1031	30	Ila	CH12022	78	Ila	CH1350	29	Ila	CH1724-003	72	I
CH1100	27	Ila	CH12023	78	Ila	CH1355	29	Ila	CH1725	74	I
CH11000	78	Ila	CH1210-1	23	Ila	CH1400	23	I	CH1726	74	I
CH11001	78	Ila	CH1211	23	Ila	CH1402	23	I	CH1726-001	74	I
CH11002	78	Ila	CH1215-1	28	Ila	CH145	15	I	CH1728	72	I
CH11003	78	Ila	CH1216	28	Ila	CH147	15	I	CH1730-001	74	Ila
CH11004	78	Ila	CH1220-1	28	Ila	CH150	15	I	CH1730-002	74	Ila
CH11005	78	Ila	CH1221	28	Ila	CH1500	61	I	CH1730-003	74	Ila
CH11006	78	Ila	CH1230	25	Ila	CH1520	66	I	CH1730-004	74	Ila
CH11007	78	Ila	CH1240-1	28	Ila	CH1560	65	I	CH1730-006	74	Ila
CH11008	78	Ila	CH1242-1	28	Ila	CH1560-MM	65	I	CH1730-008	74	Ila
CH11009	78	Ila	CH1243	24	Ila	CH1562	65	I	CH1735	73	I
CH1101	27	Ila	CH1248	23	Ila	CH1564	65	I	CH1740	72	I
CH11010	78	Ila	CH1249	23	Ila	CH1565	65	I	CH1760	72	I
CH11011	78	Ila	CH1250	23	Ila	CH1566	65	I	CH1770	73	I
CH11012	78	Ila	CH1251	24	Ila	CH1568	65	I	CH1772	74	I
CH11016	78	Ila	CH1260	24	Ila	CH1570	65	I	CH1775	75	Ila
CH11020	78	Ila	CH1261	24	Ila	CH1571	65	I	CH1780	75	I
CH11022	78	Ila	CH1262	24	Ila	CH1590	75	I	CH1800	77	I

code	page	class	code	page	class	code	page	class	code	page	class
CH1805	77	I	CH2260	39	I	CH2460	55	I	CH3391	128	I
CH1806	77	I	CH2398	47	I	CH2461	131	I	CH3392	128	I
CH1808	77	I	CH2400	47	I	CH2471	131	I	CH340	17	I
CH1809	77	I	CH2406	47	I	CH2481	49	I	CH342	17	I
CH1900	125	I	CH2406	131	I	CH2482	49	I	CH344	17	I
CH2006	45	I	CH2406-001	47	I	CH2483	49	I	CH35	14	I
CH2006-001	45	I	CH2407	47	I	CH2484	49	I	CH350	17	I
CH2006-002	44	I	CH2407	131	I	CH2492	48	I	CH36	14	I
CH2007	45	I	CH2407-001	47	I	CH2494	48	I	CH360	17	I
CH2007-001	45	I	CH2408	47	I	CH2496	48	I	CH362	17	I
CH2007-002	44	I	CH2408	131	I	CH2498	48	I	CH431	16	I
CH2008	45	I	CH2408-001	47	I	CH2500	48	I	CH440	17	I
CH2009	45	I	CH2410	55	I	CH2501	48	I	CH450	17	I
CH2010	43	I	CH2411	131	I	CH2502	48	I	CH470	16	I
CH2015	43	I	CH2415	55	I	CH2503	48	I	CH471	16	I
CH2020	45	I	CH2416	131	I	CH2504	48	I	CH5098	104	Ila
CH2020-001	45	I	CH2417	47	I	CH2506	48	I	CH5120	105	Ila
CH2025-S	44	I	CH2417	131	I	CH2507	48	I	CH5160	85	I
CH2026	39	I	CH2419	131	I	CH2508	48	I	CH5170	59	I
CH2027	39	I	CH2420	131	I	CH2509	48	I	CH5171	59	I
CH2030	43	I	CH2421	131	I	CH2510	48	I	CH5172	59	I
CH2032	43	I	CH2422	131	I	CH2520	46	I	CH5173	59	I
CH2032-S	44	I	CH2423	131	I	CH2525	46	I	CH5174	59	I
CH2034	43	I	CH2424	131	I	CH2530	46	I	CH5174-014	59	I
CH2034-001	43	I	CH2425	131	I	CH2535	46	I	CH5175	59	I
CH2035-S	39	I	CH2426	131	I	CH2600	31	I	CH5176	59	I
ch2036-s	39	I	CH2426-001	131	I	CH2605	31	I	CH5177	59	I
CH2040	40	I	CH2427	131	I	CH2610	31	I	CH5178	59	I
CH2052	39	I	CH2428	131	I	CH2611	31	I	CH5179	59	I
CH2053	39	I	CH2429	131	I	CH2612	31	I	CH5180	59	I
CH2060	40	I	CH2432	46	I	CH2615	31	I	CH5181	59	I
CH2070	40	I	CH2433	131	I	CH2616	31	I	CH5182	59	I
CH2075-S	44	I	CH2435	131	I	CH2617	31	I	CH5183	59	I
CH2080	43	I	CH2437	131	I	CH2618	31	I	CH5184	59	I
CH2085	43	I	CH2438	131	I	CH2800	128	I	CH5184-014	59	I
CH2086	45	I	CH2440	55	I	CH2802	128	I	CH5195	59	I
CH2086-S	44	I	CH2442	131	I	CH2804	128	I	CH5197	73	Ila
CH2088	45	I	CH2443	48	I	CH2806	128	I	CH5199	73	I
CH2090	43	I	CH2444	48	I	CH290	16	I	CH5200	63	I
CH2095	43	I	CH2445	48	I	CH30	14	I	CH560	18	I
CH2096	45	I	CH2446	48	I	CH31	14	I	CH5625	126	I
CH2096-S	44	I	CH2448	49	I	CH315	16	I	CH564	18	I
CH2098	45	I	CH2452	49	I	CH330	16	I	CH5650	46	I
CH2240	39	I	CH2454	49	I	CH3300	130	Ila	CH5651	35	I
CH2250	39	I	CH2456	49	I	CH3302	130	Ila	CH5651-001	35	I
CH2255	39	I	CH2457	49	I	CH3390	128	I	CH5651-002	35	I

code	page	class	code	page	class	code	page	class	code	page	class
CH5651-003	35	I	CH5674-001	38	I	CH5885	62	I	CH5954	64	I
CH5652	36	I	CH5674-002	38	I	CH5885-T	62	I	CH6000	65	I
CH5653	38	I	CH5675	38	I	CH5887	62	I	CH6001	65	I
CH5655	40	I	CH5676	38	I	CH5887-T	62	I	CH6002	65	I
CH5656	37	I	CH5677	38	I	CH5893	62	I	CH6020	65	I
CH5657	37	I	CH5677-9	39	I	CH5895	62	I	CH6021	65	I
CH5660	36	I	CH5678	39	I	CH5897	62	I	CH6022	65	I
CH5660-S	37	I	CH5679	39	I	CH5899	62	I	CH6030	80	Ila
CH5661-001	32	I	CH5680	39	I	CH590	18	I	CH6031	80	Ila
CH5661-002	32	I	CH5681	39	I	CH5900	62	I	CH6032	80	Ila
CH5661-004	32	I	CH5682	41	I	CH5900-1	62	I	CH6033	80	Ila
CH5661-005	32	I	CH5684	41	I	CH5900-T	62	I	CH6035	81	Ila
CH5662	36	I	CH5686	41	I	CH5902	62	I	CH6036	81	Ila
CH5662-S	37	I	CH5687	41	I	CH5902-1	62	I	CH6037	81	Ila
CH5663-001	32	I	CH5688	40	I	CH5902-T	62	I	CH6038	81	Ila
CH5663-002	32	I	CH5690	40	I	CH5903	61	I	CH6040	81	Ila
CH5663-004	32	I	CH5692	41	I	CH5904	62	I	CH6041	81	Ila
CH5663-005	32	I	CH5692-001	41	I	CH5904-1	62	I	CH6042	81	Ila
CH5665	36	I	CH5693	41	I	CH5904-T	62	I	CH6043	81	Ila
CH5665-S	37	I	CH5695	38	I	CH5905	61	I	CH6045	81	Ila
CH5666	36	I	CH5697	41	I	CH5906	63	I	CH6046	81	Ila
CH5666-001	36	I	CH5698-S	41	I	CH5906-T	62	I	CH6047	81	Ila
CH5666-002	36	I	CH5700	37	I	CH5907	61	I	CH6048	81	Ila
CH5666-003	36	I	CH5704	37	I	CH5908	63	I	CH6050	81	Ila
CH5666-01S	35	I	CH5705	37	I	CH5908-1	63	I	CH6052	81	Ila
CH5666-02S	35	I	CH5705-S	35	I	CH5910	63	I	CH6054	65	I
CH5666-03S	35	I	CH5706	33	I	CH5910-1	63	I	CH6055	65	I
CH5666-S	35	I	CH5708	34	I	CH5912	63	I	CH6056	65	I
CH5667-001	37	I	CH5720	34	I	CH5912-1	63	I	CH6060	81	Ila
CH5667-002	37	I	CH5720-001	35	I	CH5912-T	62	I	CH6062	81	Ila
CH5667-003	37	I	CH5725	34	I	CH5914-T	62	I	CH6065	83	Ila
CH5668-002	32	I	CH5725-001	35	I	CH5916-T	62	I	CH6066	83	Ila
CH5668-003	32	I	CH5730	34	I	CH5918-T	62	I	CH6067	83	Ila
CH5669-001	37	I	CH5730-001	35	I	CH5920	63	I	CH6068	83	Ila
CH5669-002	37	I	CH5735	34	I	CH5922	63	I	CH6070	84	Ila
CH5669-003	37	I	CH5735-001	35	I	CH5924	63	I	CH6071	84	Ila
CH5670	36	I	CH5740	34	I	CH5926	63	I	CH6072	84	Ila
CH5670-001	37	I	CH5740-001	35	I	CH5930	64	I	CH6073	84	Ila
CH5671-001	36	I	CH5741	46	I	CH5931	64	I	CH6075	84	Ila
CH5671-002	36	I	CH5742	46	I	CH5932	64	I	CH6076	80	Ila
CH5671-003	36	I	CH5750	20	Ila	CH5934	64	I	CH6077-001	80	Ila
CH5671-004	36	I	CH5751	20	Ila	CH5935	64	I	CH6077-002	80	Ila
CH5671-005	36	I	CH5752	20	Ila	CH5936	64	I	CH6078	80	Ila
CH5672	38	I	CH584	18	I	CH5940	64	I	CH6078-001	80	Ila
CH5673	38	I	CH5883	62	I	CH5950	64	I	CH6079	80	Ila
CH5674	38	I	CH5883-T	62	I	CH5952	64	I	CH6079-001	80	Ila

code	page	class	code	page	class	code	page	class	code	page	class
CH6081	84	Ila	CH6152	87	Ila	CH6444	100	Ila	CH6800	106	Ila
CH6082	84	Ila	CH6153	87	Ila	CH6446	100	Ila	CH6801	106	Ila
CH6083	84	Ila	CH6170	87	Ila	CH6461	120	Ila	CH6802	106	Ila
CH6085	84	Ila	CH6180	89	Ila	CH6462	120	Ila	CH6803	106	Ila
CH6086	84	Ila	CH6181	89	Ila	CH6464	120	Ila	CH6805	22	I
CH6087	84	Ila	CH6182	90	Ila	CH6473	121	Ila	CH6807	110	Ila
CH6088	84	Ila	CH6183	90	Ila	CH6474	121	Ila	CH6808	110	Ila
CH6089	82	Ila	CH6185	90	Ila	CH6500	92	Ila	CH6809	110	Ila
CH6090	82	Ila	CH6190	93	Ila	CH6502	92	Ila	CH6814	103	Ila
CH6091	82	Ila	CH6198	93	Ila	CH6510	91	Ila	CH6816	103	Ila
CH6092	82	Ila	CH6199	93	Ila	CH6512	91	Ila	CH6818	103	Ila
CH6093-001	82	Ila	CH6200	93	Ila	CH6516	89	Ila	CH691	18	I
CH6093-002	82	Ila	CH6201	93	Ila	CH6517	89	Ila	CH6950-001	88	Ila
CH6094	82	Ila	CH6210	108	Ila	CH6518	89	Ila	CH6950-002	89	Ila
CH6095	85	Ila	CH6211	108	Ila	CH6519	89	Ila	CH6950-003	87	Ila
CH6096	85	Ila	CH6212	108	Ila	CH6520	91	Ila	CH6950-004	88	Ila
CH6097	86	Ila	CH6215	95	Ila	CH6521	89	Ila	CH6950-005	88	Ila
CH6098	86	Ila	CH6220	113	Ila	CH6522	119	Ila	CH6950-006	88	Ila
CH6099	86	Ila	CH6225	110	Ila	CH6524	119	Ila	CH6950-007	25	Ila
CH6100	82	Ila	CH6226	110	Ila	CH6540	118	Ila	CH6950-008	42	I
CH6101	82	Ila	CH6230	110	Ila	CH6542	118	Ila	CH6950-009	64	I
CH6102	82	Ila	CH6233	103	Ila	CH6544	118	Ila	CH7000	69	I
CH6110	83	Ila	CH6234	103	Ila	CH6546	118	Ila	CH7002	69	I
CH6111	83	Ila	CH6235	103	Ila	CH6548	118	Ila	CH7100	94	Ila
CH6113	83	Ila	CH6236	114	Ila	CH6567	90	Ila	CH7102	94	Ila
CH6115	85	Ila	CH6237	114	Ila	CH6568	90	Ila	CH7110	105	Ila
CH6116	85	Ila	CH6260	99	Ila	CH6570	91	Ila	CH7112	105	Ila
CH6117	85	Ila	CH6262	99	Ila	CH6572	91	Ila	CH7115	118	Ila
CH6118	85	Ila	CH6270	99	Ila	CH6574	91	Ila	CH7117	95	Ila
CH6120	85	Ila	CH6272	99	Ila	CH6600	111	Ila	CH7118	120	Ila
CH6121	85	Ila	CH6280	100	Ila	CH6602	111	Ila	CH7119	120	Ila
CH6122	85	Ila	CH6300	105	Ila	CH6620	109	Ila	CH7120	86	Ila
CH6123	85	Ila	CH6302	105	Ila	CH6622	109	Ila	CH7122	86	Ila
CH6124	85	Ila	CH6310	105	Ila	CH6640	114	Ila	CH7124	86	Ila
CH6126	85	Ila	CH6312	105	Ila	CH6642	114	Ila	CH7140	94	Ila
CH6128	85	Ila	CH6320	93	Ila	CH6646	115	Ila	CH7155	94	Ila
CH6129	85	Ila	CH6325	93	Ila	CH6648	112	Ila	CH7170	95	Ila
CH6130	85	Ila	CH6340	119	Ila	CH6650	112	Ila	CH7175	97	Ila
CH6132	85	Ila	CH6350	119	Ila	CH6651	112	Ila	CH7180	106	Ila
CH6140	85	Ila	CH6352	119	Ila	CH6652	112	Ila	CH7182	106	Ila
CH6143	80	Ila	CH6355	115	Ila	CH6653	112	Ila	CH7184	106	Ila
CH6144	80	Ila	CH6360	105	Ila	CH6654	112	Ila	CH7186	106	Ila
CH6145	80	Ila	CH6400	121	Ila	CH6660	109	Ila	CH7188	118	Ila
CH6148	80	Ila	CH6405	121	Ila	CH6661	109	Ila	CH7190	94	Ila
CH6150	87	Ila	CH6422	112	Ila	CH6662	109	Ila	CH7192	94	Ila
CH6151	87	Ila	CH6442	100	Ila	CH680	18	I	CH7194	94	Ila

code	page	class	code	page	class	code	page	class	code	page	class
CH7202	96	Ila	CH7450	116	Ila	CH8084	127	I	CH8270	22	I
CH7204	96	Ila	CH7452	116	Ila	CH8085	127	I	CH830	19	I
CH7206	96	Ila	CH7454	116	Ila	CH810	19	I	CH832	19	I
CH7210	97	Ila	CH7456	99	Ila	CH8105-001	21	I	CH8320	126	I
CH7212	97	Ila	CH7458	116	Ila	CH8105-002	21	I	CH8360	34	I
CH7214	97	Ila	CH7460	117	Ila	CH8105-003	21	I	CH840	19	I
CH7216	97	Ila	CH7470	117	Ila	CH8105-004	21	I	CH841	19	I
CH7220	95	Ila	CH7472	117	Ila	CH8105-005	21	I	CH8505	83	Ila
CH7221	95	Ila	CH7474	117	Ila	CH8105-006	21	I	CH8510	83	Ila
CH7222	95	Ila	CH7480	117	Ila	CH8110	21	I	CH8515	58	I
CH7230	115	Ila	CH7490	99	Ila	CH8112	21	I	CH8520	58	I
CH7231	116	Ila	CH7492	99	Ila	CH8114	21	I	CH8525	58	I
CH7232	116	Ila	CH7510	98	Ila	CH8119	20	I	CH8530	58	I
CH7240	115	Ila	CH7512	98	Ila	CH8120	20	I	CH8535	58	I
CH7245	98	Ila	CH7600	111	Ila	CH8121	20	I	CH8540	58	I
CH7247	98	Ila	CH7602	111	Ila	CH8123	20	I	CH8545	32	I
CH7251	102	Ila	CH7813	107	Ila	CH8124	20	I	CH8550	32	I
CH7252	102	Ila	CH7814	107	Ila	CH8125	20	I	CH8555	33	I
CH7262	107	Ila	CH7815	107	Ila	CH8126	20	I	CH8560	33	I
CH7264	107	Ila	CH7902	92	Ila	CH8127	20	I	CH8580	128	I
CH7265	108	Ila	CH7910	92	Ila	CH8130	21	I	CH8585	128	I
CH7268	108	Ila	CH7912	93	Ila	CH8131	21	I	CH8586	127	I
CH7270	108	Ila	CH7950	122	Ila	CH8152	20	I	CH8589	49	I
CH7272	107	Ila	CH7952	122	Ila	CH8153	20	I	CH8590	50	I
CH7274	107	Ila	CH7970	121	Ila	CH8154	20	I	CH8590-001	50	I
CH7276	107	Ila	CH7972	121	Ila	CH8155	20	I	CH8591	49	I
CH7280	99	Ila	CH7974	121	Ila	CH8156	20	I	CH8591-001	50	I
CH7284	96	Ila	CH7974-001	121	Ila	CH8158	22	I	CH8592	50	I
CH7286	96	Ila	CH7975	121	Ila	CH8180	55	I	CH8592-001	50	I
CH7288	96	Ila	CH80	14	I	CH8181	55	I	CH860	24	I
CH7290	113	Ila	CH800	19	I	CH8182	55	I	CH8608	69	I
CH7292	113	Ila	CH8040	129	Ila	CH8183	55	I	CH8610	69	I
CH7294	113	Ila	CH8042	129	Ila	CH8184	55	I	CH8620	69	I
CH7296	113	Ila	CH8044	130	Ila	CH8185	53	I	CH8622	69	I
CH7298	113	Ila	CH8046	130	Ila	CH8186	53	I	CH8669	22	Ila
CH7300	104	Ila	CH8047	130	Ila	CH8187	53	I	CH8670	22	Ila
CH7302	104	Ila	CH8048	130	Ila	CH8188	53	I	CH8672	125	I
CH7304	104	Ila	CH8049	130	Ila	CH82	14	I	Ch8673	125	I
CH7310	120	Ila	CH805	19	I	CH8210	53	I	CH8675	43	I
CH7340	116	Ila	CH8050	129	Ila	CH8212	53	I	CH8676	31	I
CH7350	101	Ila	CH8055	129	Ila	CH8214	53	I	CH8680	126	I
CH7351	101	Ila	CH8057	129	Ila	CH8216	53	I	CH8680-001	126	I
CH7355-001	102	Ila	CH8059	129	Ila	CH8218	53	I	CH8681	127	I
CH7355-002	102	Ila	CH8060	129	Ila	CH8220	53	I	CH8689	42	I
CH7355-003	102	Ila	CH8080	128	Ilb	CH8225	53	I	CH8690	42	I
CH7355-004	102	Ila	CH8083	128	Ilb	CH8230	53	I	CH8752	53	I

code	page	class
CH8756	53	I
CH8766	53	I
CH8770	67	I
CH8770-001	67	I
CH8772	67	I
CH8800	49	I
CH900	19	I
CH90-1	15	I
CH905	19	I
CH9300	14	I
CH9302	14	I
CH9304	14	I
CH9306	14	I
CH9308	14	I
CH9312	14	I
CH9314	14	I
CH9316	14	I
CH9318	14	I
CH9320	14	I
CH9322	14	I
CV1033	123	Ila
CV1061	123	Ila
CV1133	124	Ila
CV1161	124	Ila
CV1186	124	Ila
CV5020	122	Ila
CV5025	122	Ila
CV5030	122	Ila
CV5040	122	Ila
CV5045	123	Ila
CV5050	123	Ila
CV5055	123	Ila
CV5110	122	Ila
CV5111	122	Ila

carefusion.com



CareFusion
1500 Waukegan Road
McGaw Park, IL 60085 USA



CareFusion France 309 S.A.S.
8 bis rue de la Renaissance
F-44110 CHATEAUBRIANT

Please read the complete Instructions for Use that come with the instruments or follow the instructions on the product labeling.

Medical devices class I, I measuring, IIa and IIb according to Medical Devices Directive 93/42/EC - CE-marked by TÜV Product Service GmbH (0123) where appropriate.

©2012 CareFusion Corporation or one of its subsidiaries. All rights reserved.

Delicate Touch, Diamond Jaw, Genesis, Mueller Mirror, Snowden-Pencer, V. Mueller and Vital are trademarks or registered trademarks of CareFusion Corporation or one of its subsidiaries.

LIT12012en-ED01-Dec2012



CareFusion